syngo.share / Release VA30C / 2021-12-16 / Revision 6272

syngo.share view

Online Help

Language Directory

Deutsch	3
English	96
Français	180

syngo.share / Release VA30C / 2021-12-16 / Revision 6272

syngo.share view

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Online Help (English)

Contents

1	Introduction	101
1.1	Intended Use	101
1.2	Technical Terms and Notation	101
1.2.1	Generic Document	101
1.2.2	DICOM Document	101
1.2.3	DICOM Information	101
1.2.4	EXIF Information	102
1.2.5	Modality	102
1.2.6	Organizational Unit	102
1.2.7	Hounsfield Unit	102
1.3	Safety Notices	103
	Image and Result Distribution	103
	Medical Education	103
	Accuracy of Calculations	103
	Units	103
	Access Restriction	103
	Effects of External Programs	104
	Printing of Images	104
	Network Connection	104
	Workstation Accessibility	104
2	Settings	105
2.1	Saving and Loading a Profile	105
2.2	General	106
2.2.1	Program-specific Options	106
2.2.2	Enabled User Interfaces	106
2.2.3	Data Sources	106
2.2.4	Viewer	106
2.3	Display	108
2.3.1	Viewer	108
2.3.2	Results List	109
2.3.3	Loaded Documents	109
2.3.4	Annotations	109
2.3.5	Presentation States	110
2.3.6	Drag and Drop	110
2.4	Hanging	110
2.5	Labels	113
2.6	Communication	114
2.6.1	<i>syngo.share core</i>	114
2.6.2	DICOM SCU	114
2.6.3	DICOM SCP	114
2.6.4	XDS Consumer Adaptor	115
2.7	Advanced	115
2.7.1	Memory Management	115
2.7.2	Separators for Windowing Default Settings	115
2.7.3	Modality-Depending Defaults	115
2.7.3.1	Additional information	116
2.7.4	External Applications	116
2.7.5	Export Application	117
2.7.6	icoscript	118

2.8	Print	118
2.8.1	System Printer	118
2.8.2	DICOM Printer	119
2.8.3	DICOM Printer Settings	119
2.8.4	Print Layout	120
2.9	DICOM Export Profiles	121
3	Usage Instructions	122
3.1	Getting Started	122
3.1.1	User Interface	122
3.1.1.1	Main Menu	122
3.1.1.2	Main Area	123
3.1.1.3	Document Tree	123
3.1.2	Display Modes	123
3.1.3	Login and Logout	123
3.1.4	Changing the Password	124
3.2	Searching and Loading	124
3.2.1	Searching in Archives	125
3.2.1.1	Configuring the Search Fields	125
3.2.1.2	Types of Search Fields	126
3.2.1.3	Searching by Keyword	127
3.2.1.4	Saving a Search Query	128
3.2.1.5	Results List	129
3.2.1.6	Using the Series Area in the Archive/Entity/Collections Search	131
3.2.1.7	Using the Series Area in the DICOM Q/R and External Search	133
3.2.2	Searching Local Storage Media	133
3.2.3	Loading Search Results	134
3.2.4	Document Tree	134
3.2.4.1	Basic Actions	135
3.2.4.2	Document Tree Main Menu	135
3.2.4.3	Document Tree Context Menu	136
3.2.4.4	Keyword Tagging	137
3.2.4.5	Collections	138
3.3	Image Display and Editing	139
3.3.1	Hanging Documents	139
3.3.1.1	Hang	139
3.3.1.2	Hang 2x2	140
3.3.1.3	Series Overview	140
3.3.1.4	Compare Mode	140
3.3.1.5	Hang Structured Report	142
3.3.1.6	Action before Loading	142
3.3.1.7	Action after Loading	142
3.3.1.8	Applying a Hanging Workflow	142
3.3.2	The Viewer	142
3.3.2.1	Dividing the Viewer Vertically and/or Horizontally	143
3.3.2.2	Defining the Layout of Individual Viewers	143
3.3.2.3	Viewer Context Menu	144
3.3.2.4	Toolbar	144
3.3.2.5	Toolbox	144
3.3.2.6	Tool Region	145
3.3.2.7	Presentation States (PR)	145
3.3.3	General Functions	146

3.3.3.1	Assign Properties	146
3.3.3.2	Calculations, Text Annotations, Drawing	146
3.3.3.3	Camera	150
3.3.3.4	Clear Viewer	151
3.3.3.5	Density	151
3.3.3.6	Image Administration	151
3.3.3.7	Integration	151
3.3.3.8	Invert	152
3.3.3.9	Layout	152
3.3.3.10	Mode Selection	152
3.3.3.11	Pan	152
3.3.3.12	Playback Controls	152
3.3.3.13	Pointer	153
3.3.3.14	Print	153
3.3.3.15	Pseudo Colors	153
3.3.3.16	Rotate and Flip	153
3.3.3.17	Scroll	154
3.3.3.18	Series Synchronization	154
3.3.3.19	Sorting Criteria	154
3.3.3.20	Switch Screens	154
3.3.3.21	Tag	155
3.3.3.22	Viewing Properties	155
3.3.3.23	Undo Editing	156
3.3.3.24	Level/Window (DICOM only)	156
3.3.3.25	Zoom	157
3.3.4	ECG Mode	158
3.3.5	Movie Mode	159
3.3.6	Video Mode	160
3.3.7	Metadata Mode	160
3.3.7.1	Searching Metadata	161
3.3.7.2	Copying Metadata	161
3.4	Document Management	161
3.4.1	Merging Patients	161
3.4.1.1	Finding Created Patients	161
3.4.1.2	Finding a Matching Patient in the Patient Management System	162
3.4.1.3	Determining the Correct Patient	162
3.4.1.4	Saving or Discarding Changes	162
3.4.2	Changing the Description	163
3.4.2.1	Study Level: Changing the Document Properties	163
3.4.2.2	Series Level: Changing the Series Description	163
3.4.2.3	Image Level: Changing the Image Description	163
3.4.3	Printing Documents	163
3.4.4	Printing Images	164
3.4.4.1	Printing Directly	164
3.4.4.2	Printing via the Print Layout	164
3.4.5	Exporting Documents	166
3.4.5.1	Exporting Documents as DICOMDIR	166
3.4.5.2	Exporting Documents to a Storage Medium	168
3.4.5.3	Errors	170
3.4.5.4	Creating a Dossier	170
3.4.6	Anonymizing Documents	171
3.4.6.1	Anonymizing Documents during Export	171

3.4.6.2	Anonymizing Documents during Printing	171
3.4.7	Moving Documents	172
3.4.8	Copying Documents	172
3.4.9	Sending Documents	173
3.4.10	Sharing Documents	174
3.4.11	Registering Reference Pointers	175
3.4.12	Reassigning Documents	175
3.4.13	Deleting Documents	176
3.4.14	Viewing and Correcting Warnings and Errors	177
3.4.15	List of Favorites	177
3.5	Keyboard Shortcuts <i>syngo.share view</i>	178

1 Introduction

1.1 Intended Use

syngo.share view is software to be used in clinical image and result distribution and provides access to clinical data from connected systems (such as Vendor Neutral Archives, Picture Archiving and Communication Systems, DICOM archives, XDS-I Affinity Domains, Electronic Medical Record Systems, Clinical Information Systems) and from local filesystems for the electronic display of data. This product is not intended for diagnosis, prevention, monitoring, prediction, prognosis, treatment, cure or alleviation of diseases or other medical conditions.

1.2 Technical Terms and Notation

1.2.1 Generic Document

This term refers to any non-DICOM file that can be viewed in or archived by *syngo.share*. Non-DICOM files include images produced by digital cameras or imaging applications (e.g. JPEG, MPM, TIFF) as well as various text, video, and audio file formats (e.g. PDF, DOC, MPEG4, MP3) produced by various external applications (e.g. Microsoft Office, Adobe Acrobat).

Generic files are wrapped in generic containers, analogous to DICOM studies.

1.2.2 DICOM Document

A DICOM document consists of at least one DICOM study, series and instance.

DICOM documents are produced by medical imaging devices compliant with the DICOM standard. They can be viewed in or archived by *syngo.share*.

Study

A DICOM study contains all the data of a medical examination during which image data is acquired. This includes the actual image data (e.g. X-ray, CT, ultrasound) as well as the patient data and descriptions relevant to the examination (such as the examination date or type). A DICOM study contains 1-n DICOM series.

Series

A DICOM series usually equates to a specific type of data acquired/produced by a particular device or to the position of the patient during acquisition. A series contains 1-n DICOM instances.

Instance

A DICOM object instance (or simply DICOM instance) is the actual data captured by the device, such as an image or waveform object.

1.2.3 DICOM Information

DICOM (Digital Imaging and Communications in Medicine) is a standardized file format used in medical imaging.

In addition to object data (typically image data), the DICOM format contains additional information such as patient name, admission date, device parameters, and physician's name. The DICOM format also contains the information necessary for data exchange between DICOM-capable devices.

1.2.4 EXIF Information

EXIF (Exchangeable Image File Format) is a standardized data format generated by, e.g., digital cameras when taking a picture.

In addition to the image data commonly contained in JPEG or TIFF format, the EXIF format contains additional information describing the image in more detail, such as date and time information or the type of device used for image capture.

1.2.5 Modality

An imaging modality is a device for acquiring medical images such as X-rays or CTs. In *syngo.share* view, modalities are identified by abbreviations such as CR for Computed Radiography (e.g. when searching).

1.2.6 Organizational Unit

An *organizational unit* is an association of arbitrary actors in a medical organization, e.g. a department, medical devices, or groups of people. Organizational units can be customized and expanded with more constituents as needed, and these constituents can be uniformly managed.

All data to be archived is allocated to exactly one organizational unit.

1.2.7 Hounsfield Unit

Density calculations of body tissue with the density calculation tool use the Hounsfield scale. The following equation is used to calculate the attenuation coefficients of various body tissues in CTs. It is assumed that the absorption value of tissues and the physical density have a linear dependence (μ = linear attenuation coefficient):

$$HU = \left(\frac{\mu_{Object} - \mu_{Water}}{\mu_{Water}} \right) \times 1000$$

Table 1: Relative Attenuation Coefficients of Body Tissues and Water

Coefficient	Object
0	Water
1000 bis 100	Bone
100 bis 0	Liver/Muscle/Kidney
-100 through -200	Fat
-1000	Air

The relative attenuation values can be determined easily from any point or area of the image. Although it is technically possible to render any number of attenuation values in grayscale, the human eye is not equipped to distinguish between more than 20 shades of gray. Instead, the use of windowing is recommended in order to visually highlight different kinds of body tissue. ³

Please read the safety notices, see [Section 1.3](#).

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1.3 Safety Notices

Image and Result Distribution

syngo.share view is only intended for clinical image and result distribution.
Do not use *syngo.share* view for direct diagnoses.

Medical Education

The use of *syngo.share* view should be granted to qualified staff members only.

Accuracy of Calculations


Images may contain calibrations. Please note that the accuracy of these calibrations correlate to the precision of the program or modality with which they were created. Calibrations must be checked according to the manufacturer's instructions. Calculations, calibrations and results derived therefrom must be checked for plausibility by clinical professionals. Inaccuracies may arise during image capture by clinical professionals or due to the image capturing technology of a given modality. For instance, X-ray images may contain inaccuracies due to the fact that the x-rayed body part is projected onto a two-dimensional imaging surface. Inaccuracies during visualization may be due to hardware settings.

This product rounds calculations as follows: Length and area: if pixel spacing ≥ 0.1 mm: 1 decimal place, if pixel spacing < 0.1 mm: 2 decimal places; density (HU): no decimal place; density average (HU): 1 decimal place; density area (HU): no decimal place. This product uses an interpolation algorithm for image display. This product may not be able to calculate the area of polygons whose edges overlap; this will be indicated by the value "n/a".


Units

The application *syngo.share* view conforms to Directive 80/181/EEC, which contains specifications related to units of calculation. Calibrated lengths are calculated using millimeters. Surface areas are calculated in square millimeters (mm^2). Plane angles are calculated in degrees. Densities are calculated using the industry-standard unit, Hounsfield units (HU), which are not regulated by Directive 80/181/EEC.

Access Restriction

Be sure to log out of *syngo.share* view when finished working in order to avoid unauthorized access to the system.

Effects of External Programs

The software manufacturer does not guarantee that third-party applications will not negatively affect the system. The standards of any external programs accessed by *syngo.share* view are set by their respective manufacturers.

Printing of Images

When printing images using the image caption option "As camera", it is possible that some labels will be hidden from view due to the physically limited space on the page. As a result, patient information may not be visible on the printout. Possible transformations such as mirroring, rotation, and applied presentation states are not marked explicitly in the printout. If a series has active presentation states in both the view and print layouts, deactivation of the presentation states in the view layout will also deactivate it in the print layout of the series.

Network Connection

syngo.share view requires an active network connection. For IT security reasons, a (temporary) failure of the network connection may result in you being logged out and having to log in again.


Workstation Accessibility

When using *syngo.share* view, make sure that you can operate the application at all times.

2 Settings


You can find the configuration menu by clicking on  in the main menu.


Note

It is possible that the settings are not visible for you at all or that some sections in the settings are invisible or grayed out. This means that you are currently not authorized to change them. Please contact your system administrator if you still want to make changes.

You can either browse the settings by clicking on the individual categories or use the search function to quickly find particular entries. You can reset/undo your search by clicking the  in the search field.


A star next to a list category means that the settings in this section have changed. To save your changes, click .

You may have to restart or log out of/log in to the application again in order to implement some settings. In these cases, you can either agree to an immediate restart or continue working with the previous settings (the changes become effective at the next program start).

Note



Changes apply to the workstation where they are made. Therefore, never change settings randomly, especially if the same workstation is used by multiple persons.


Once you confirm a change using the  button, the change is saved locally in the Windows registry and on the server for this workstation. When the program (re)starts, the settings from the server are automatically applied. If you notice any alterations, it means that the system administrator has changed the settings for this workstation on the server side.

2.1 Saving and Loading a Profile

In the configuration menu, you can save and load your settings as a personal user profile.

This is how to create your personal profile:

1. Click  in the main menu.
2. Change the settings for your profile.
3. Click the  button.
4. Create a file name and target folder for your profile in the **Save profile** window and click . All settings are saved to a profile file with the .ics (icoserve Config File) ending.

Note

Save the profile on an external media device or network drive if you would like to have access from another workstation.

Follow these steps to load a saved profile:

1. Click  in the main menu.
2. Click the  button.

3. Choose the directory and file for the desired profile in the **Load profile** window and click `Open`. *syngo.share* view adopts all settings of the respective profile.

2.2 General

2.2.1 Program-specific Options

Language

syngo.share import is available in different languages. Select the desired language for the user interface for the logged-in user here. The entry `system language` applies the language settings of the Windows operating system.

Maximize on Startup

Enabling this option will automatically maximize the application to fullscreen when starting the program.

Multiscreen Mode on Startup

Enabling this option causes *syngo.share* view to immediately switch to multiscreen mode when starting the program.

Minimize to Taskbar

If you activate this option, the application will close as usual, but will remain in sleep mode in the background. The application icon is visible in the Windows taskbar. The sleep mode allows for a quicker startup of the application.

Automatic Login

If this option is activated by the administrator in *syngo.share* Configuration, each user who has logged on to the Windows operating system is automatically logged in to the application. If this setting is not available, you can still log in manually with your user credentials. It is also possible to use Kerberos for automatic login.

Theme

Here, you can set the color scheme for *syngo.share* view.

2.2.2 Enabled User Interfaces

You can activate/deactivate certain parts of the user interface. Upon deactivation, the corresponding icon disappears from the main menu. You can also (un)hide the patient QR code with `Show QR Code controls`.

2.2.3 Data Sources

Here you can customize which data sources should be available for searching. You can also select whether certain local directories should be exempted from the filesystem search.

2.2.4 Viewer

Interpolation Algorithm

This provides a selection list where you can choose the way the images are calculated for display.

↳None

When zooming in, the nearest neighbor method is used for display. Brightness values are not interpolated, but neighboring pixels are copied.

↳Bresenham

Default setting; benefit: faster calculation.

↳Bilinear

Benefit: optimum image display, drawback: slower than Bresenham.

↳Smoothscale

Benefit: optimum image display if the zoom factor of an image is below 100%.

☰ Tool selection

You can select viewing tools from the context menu or the toolbox by right-clicking in the Viewer. Choose the desired type of menu from the selection list. See [Section 3.3.2.3](#), [Section 3.3.2.5](#).

☰ Standard tool

Define which tool is active in the Viewer by default.

☰ View multiframe series in

This provides a selection list where you can choose whether multiframe series should be shown in movie mode or stack view when dragged into the Viewer. The default setting is movie mode.

☰ Show Viewer tools

Show or hide the Viewer toolbar.

☰ Show global tools

Show or hide the general toolbar in the Viewer area.

☰ Toolbar icon size

Adjust the size of the toolbar icons (large or small).

☰ Autoplay multiframe images

Selecting this checkbox will immediately play multiframe sequences when they are dragged from the loaded documents area to the Viewer.

☰ Always use mouse wheel for scrolling

Selecting this checkbox will allow you to use the scroll function of the mouse wheel regardless of the tool selected.

☰ Show text annotations

Selecting this checkbox will automatically show text contained in DICOM images. For editing individual images, you can enable or disable this function from the context menu of the Viewer

Viewing properties > Show text.

☰ Show DICOM overlays

Enabling this option will automatically enable the show overlays viewing tool that shows overlay pixel data when hanging DICOM images. See also [Viewing Properties](#).

☰ Activate localizer automatically

Enabling this option will automatically enable the show localizer viewing tool when layered series are hung. See also [Viewing Properties](#).


Always interpolate

Activating this option is only relevant when switching the interpolation algorithm to Bilinear or Smoothscale. This causes images to be recalculated using the selected interpolation algorithm during every editing step. It keeps the image quality high, but requires more processing power. Disabling the option here in the configuration menu (globally) or in the context menu of the Viewer (for an individual image) causes images to be calculated using the faster Bresenham interpolation algorithm when you use a viewing tool. The image is not calculated using the selected interpolation algorithm until you complete an image editing step (e.g., while you scroll the images of a series and at the end when you stop scrolling). You can also use   from the context menu of the Viewer to control this.

SOPClassUIDs for automatic view properties synchronization

For the SOPClassUIDs entered, the function **Synchronize Viewer** settings is automatically enabled during hanging. The following SOPClassUIDs are entered by default.

- CT (normal and advanced)
- MR (normal and advanced)
- PT
- NM

You can also enable or disable this function via the context menu under   for editing individual images.

Use Upper/Lower Threshold for

Define image types for which the Upper/Lower threshold windowing should be used for (e.g. NM, PT). By default the Level-/Width windowing is used, see [Section 3.3.3.24](#).

Limit window to min/max

The DICOM standard provides for a grayscale range of -65500 to +65500 for Level/ Width control. However, many DICOM images do not utilize this maximum range because their brightness values are limited to a relatively small range. For these types of images, selecting this option will set limits which restrict the Level/Width control values to the brightness values of the respective image. As a result, the windowing controls display a smaller window area and are easier to manipulate. Therefore, this option is enabled by default. [Section 3.3.3.24](#) contains some additional information about this topic.

PDF Viewing, CDA Viewing, Video Viewing

Use these separate configuration options to select whether the associated files should be opened in *syngo.share* view (internally) or using third-party software (externally). If you choose external viewing, the files in question will be opened with the standard program configured on the local machine.

2.3 Display

Here you can adjust the settings for the visual presentation of various elements (Document Tree, calculations, etc.) in *syngo.share* view.

2.3.1 Viewer

Main screen

This is the primary viewing monitor in multiscreen mode. This setting is particularly useful when working with multiple screens of different sizes. You can specify the number of the main screen here.

☰ Text size

Choose your desired text size from the selection list of preset sizes. The entry of intermediate sizes is not supported. The new setting takes effect as soon as you reopen a series in the Viewer. Alternatively, you can use **Ctrl+** to increase and **Ctrl-** to decrease the text size in the Viewer.

☰ Text color

Click the color bar below text color to open the color selection dialog box. Select the desired text color and click **OK** to confirm.

2.3.2 Results List

Here you can adjust the display settings for search results:

☰ Maximum number of search results

Set the maximum number of search results you want to see displayed.

☰ Colorize entries

This option facilitates the selection of search results from a mixed results list by alternately highlighting entries for different patients in color.

☰ Show all marked entries entirely

This option affects how matching values in a grouping are displayed. If deactivated, this option will show matching values only for the topmost entry in a chosen grouping. If activated, it will display all values of a subordinate group member when you click on it.

2.3.3 Loaded Documents

Here you can change the size of the preview images shown in the Document Tree.

2.3.4 Annotations

☰ Sensitivity

You can choose how precisely the mouse pointer has to be positioned in order to select a given annotation. If you choose the high setting, you will have to move the mouse very close to the desired annotation; at a medium or low setting, you can successfully select the annotation from a distance.

☰ Text Size

Choose the default text size for annotation text from a list of presets.

☰ Line Width

Adjust the line width of drawings by selecting from a list of presets.

☰ Annotation Color

Click on the color bar to choose the color for annotations from a color selection window.

☰ "Show annotation text" initial on

Activating this option will show the text for all annotations by default.

☰ Draw Textbox

Selecting this option will show a box around the annotation text.

2.3.5 Presentation States

☰ Line width (Pixels)

Choose the line width for annotations from Presentation States from the selection list.

☰ Text size (Points)

Select the text size for labels from Presentation States from the list of presets.

☰ Annotation color

This option determines how *syngo.share view* will display graphical annotations from Presentation States when the **Ignore Presentation State colors** option is activated.

☰ Ignore Presentation State colors

Activating this option will ignore the color for annotations from the Presentation States and instead use the **Annotation color** setting.

2.3.6 Drag and Drop

☰ Mode to Determine Image Size

In the drop-down menu, the size of the dragged image can be chosen. Depending on the selection, width and/or height can be adjusted.

☰ Background Color

Pick the background color for the dragged image.

☰ Text Color

Pick the text color for the dragged image.

☰ Show Text Overlays

This option determines whether texts are shown or hidden in the dragged image.

☰ Text Overlays Inside Image

This option determines whether text is shown within the dragged image or moved towards the margin as much as possible.

☰ Show Annotations

This option determines whether annotations in the dragged image are shown or hidden.

☰ Show Annotation Text

This option determines whether annotation text is shown or hidden in the dragged image.

2.4 Hanging

You can configure so-called Hanging Workflows in order to automatically execute certain actions when hanging documents in the Viewer. A workflow consists of one or more steps; which are in turn comprised of conditions and actions.

Use the **Definitions** tab to set conditions and actions, which can then be combined into steps:

☰ Create/Define Condition

Conditions define when a particular workflow is shown in the context menu of the Document Tree, or rather, to which of the loaded documents a workflow's component action(s) can be applied. You can search for existing conditions or create/delete your own. In order to create a new condition, click the to enter a name for it. Then set the object to which the condition should apply (e.g. the DICOM modality), as well as its parameters (e.g. only OPs/OPTs).

☰ Create/Define Action

Here, you can create actions which will be automatically executed as part of a workflow. You can either search for existing actions or create/delete your own. In order to create a new action, click , enter a name, and then choose the desired action from the drop-down menu in the **Define Action** area. Certain actions require additional parameters to be set (indicated by their validity status). For example, the "Gradient" action cannot be executed without first setting the desired color bar.

In the **Usage** tab, you can then combine actions and conditions into steps, and steps into workflows:

☰ Create/Define Steps

Steps consist of one or more actions and their set conditions. Steps are executed automatically when a workflow to which they are assigned is started. In order to create a step, click on to name it. Then set the desired actions and the conditions under which they should be taken in the **Define Step** area (e.g. DICOM-type OP/OPT images should be shown in gradient colors upon hanging).

☰ Create Workflow/Add Steps

Here, you can combine any number of steps into a workflow. Create a new workflow by clicking the and entering a name. Then add the desired predefined steps in the **Add Steps** area. Certain steps can be set as repeatable. A workflow is available in the context menu of the Document Tree (in the **Load** and **View** areas) if at least one of its conditions applies to the loaded document(s).

Example Workflow: Auto-Hanging of Multiframe Instances

Hanging workflows can be used e.g. to automatically hang all (multiframe) instances of an ultrasound series on viewports arranged in a 2x2 layout. To set up this example workflow follow these steps (new entries in each panel are created by clicking the button):

1. In **Settings** open the section **Hanging** and go to the tab **Definitions**.
2. Create a condition and enter a speaking name, e.g. "is US".
3. Define the just created condition in the next panel. As property choose "Modality", as condition "equals", and type the value "US".
4. Create four actions in the next panel, e.g. "US split top left", "US split top right", "US split bottom left", and "US split bottom right".
5. Select each action one after another in the **Create action** panel and define it in the **Define action** panel. For "US split top left" add the action "Screen Layout" and for parameter choose a 2x2 layout. Add another action "Instance viewport" and for parameter select the top left viewport. Select the next actions in the **Create action** panel and define the corresponding actions, namely setting the correct instance viewport for each. It is not necessary to set the screen layout for each action.
6. Switch to the **Usage** tab and create a step in the first panel, e.g. "hang US".
7. In the **Define step** panel add the condition "is US" four times and assign the actions created from top left to bottom right one at a time.

8. Create a new workflow in the respective panel and enter a speaking name, e.g. "hang US workflow".
9. In the **Add steps** panel select the "hang US" step created above, and check the box for "repeatable".
10. The hanging workflow is created.

The context menu of every study that fulfills at least one of the configured conditions (for the example that means at least one US series exists), now shows the entry "Perform hang US". On clicking this entry all the instances of the study will be hung to viewports into the defined 2x2 layout if one of the series of the selected study meets the condition "Modality equals US". Use the ◀ and ▶ buttons next to the workflow name to perform the next workflow step, i.e. in this example switch to the next or previous four instances.

Example Workflow: Auto-Hanging of Ophthalmic Images

Hanging workflows can be used e.g. to automatically hang ophthalmic tomography [OPT] instances below the corresponding ophthalmic photography [OP] instance of the same study. To set up this example workflow follow these steps (new entries in each panel are created by clicking the + button) (Table 2 gives a tabular overview of the steps):

1. In **Settings** open the section **Hanging** and go to the tab **Definitions**.
2. Create a condition and enter a speaking name, e.g. "is OPT".
3. Define the just created condition in the next panel. As property choose "Modality", as condition "equals", and type the value "OPT". In most cases, a second definition is needed for a clear distinction. Since only OPT series containing image data can be referenced to an OP instance, in this example, OPT series with image data have to be distinguished from OPT series with other data, e.g. DICOM-wrapped PDF files. This distinction can be achieved via the series description: Only OPT series containing image data have a certain value as part of the series description, here "IR".
Hence, for this example, add another definition, as property choose "SeriesDescription", as condition "contains", and type the value "IR".
4. Create one action in the **Create action** panel, e.g. "show referenced image at top".
5. Define the action in the next panel. Add the action "Screen Layout" and for parameter choose a 1x2 layout. Add another action "Series viewport" and for parameter select the second viewport from the top in the left column. Add another action "Show referenced image" and for parameter select the top left viewport.
6. Switch to the **Usage** tab and create a step in the first panel, e.g. "hang OP".
7. In the **Define step** panel add the condition "is OPT" and assign the action "show referenced image at top" created before.
8. Create a new workflow in the respective panel and enter a speaking name, e.g. "hang OP workflow".
9. In the **Add steps** panel select the "hang OP" step created above, and check the box for "repeatable".
10. The hanging workflow is created.

Table 2: Example for Auto-Hanging of Ophthalmic Images

Tab	Section	Entry
Definitions	Create Definition	"is OPT"
Definitions	Define Condition	Modality - equals - OPT SeriesDescription - contains - IR
Definitions	Create Action	"show referenced image at top"
Definitions	Define Action	Screen Layout - 0,0,2,1,0 - ✓ Series viewport - 1,0,0 - ✓ Show referenced image - 0,0,0 - ✓
Usage	Create steps	"hang OP"
Usage	Define step	is OPT - show referenced image at top
Usage	Create workflow	"hang OP workflow"
Usage	Add steps	hang OP - ✓

The context menu of every study that fulfills at least one of the configured conditions (for the example that means at least one OPT series exists), now shows the entry "Perform hang OP". If one of the series of the selected study meets the conditions "Modality equals OPT" and "SeriesDescription contains IR" on clicking this entry all the OPT instances of the study will be hung below the corresponding OP instance in the defined 1x2 layout. Use the  and  buttons next to the workflow name to perform the next workflow step, i.e. in this example switch to the next or previous two instances.

2.5 Labels

You can customize the display of descriptive labels in the Document Tree and/or various areas of the Viewer.

Labeling uses DICOM tags, e.g. entering the tags <SeriesDescription><Modality> in the DICOM series field in the **Loaded Documents** area will display the modality used to produce the image(s) in the series description of the Document Tree.

You can also use basic HTML tags for styling labels, e.g.
 for a line break, for text in bold face, etc.


The following labels can be configured:

- Loaded Documents
- Text Overlays in Default Viewer
- Text Overlays for MPR Viewer
- Text Overlays for ECG Viewer

Note

Please note that certain values, e.g. in the ECG Viewer, are composed of multiple DICOM tags. For instance, the PQ Interval is created by combining the CodingSchemeDesignator, CodingSchemeVersion and CodeValue tags.

Note

If you want certain information to always be displayed, the corresponding variable should be repeated in the entry fields for the small or very small label. This causes the value to be displayed regardless of the available image space, be it limited or very limited.

2.6 Communication

In this section, you can find settings for communication with the system. These settings should be made or checked by the system administrator.

2.6.1 *syngo.share* core

Enter the details for communication with *syngo.share* core here.

Host

Name of the server where the naming service of *syngo.share* core is available.

Port

Listen port of the naming service of *syngo.share* core.

Meta

Service name of the metaserver that is to be used.

Event

Service name of the event server that is to be used.

Add for DICOM Query/Retrieve

This switch adds *syngo.share* core to the DICOM Q/R tab in the **Load** area. The results list of a DICOM Q/R search is independent of the standard patient search.

Shared Calling AET

Application Entity Title used by *syngo.share* view to identify itself at the host.

Automatically Configure Move Destination

This switch activates the automatic configuration of a destination of move operations.

Listen Port Range

Port range used to listen for incoming DICOM communication.

2.6.2 DICOM SCU

Here you can set the details of *syngo.share* view as the ServiceClassUser.

Calling AET

Here, you can specify which Calling AET is used to contact remote DICOM devices.

TCP Listen Port

This is the unencrypted port used for C-MOVE operations.

TLS Listen Port

This is the encrypted port used for C-MOVE operations.

2.6.3 DICOM SCP

Here, you can configure all the settings for communication with external PACS (ServiceClassProvider). If an external PACS is defined, *syngo.share* view also shows the DICOM Q/R tab in the **Load** area. The results list of a DICOM Q/R search is independent from the results list of the standard patient search.

2.6.4 XDS Consumer Adaptor

In order to be able to start a cross-enterprise search (XDS search) for external patients and documents in the **External** search tab, you must configure an eHealth Solutions Consumer Adaptor. The **External** tab is only available if you configure the eHealth Solutions Consumer Adaptor. The default Assigning Authority of the site is used as the Assigning Authority for the local patient ID.

☰ Consumer URI

Here, you can specify the URI of the eHealth Solutions Consumer Adaptor.

2.7 Advanced

2.7.1 Memory Management

☰ Cache Directory

Determines where caching data is stored.

☰ Download Threshold

You can choose the number of images to be loaded initially for every series of a selection list. When this threshold is reached, the load operation is interrupted. This reduces the waiting time when loading studies with multiple series from the archive – you can immediately remove unneeded series and load only the series that are needed. To reload the remaining images of a series, click in the Document Tree. No threshold is used as the default setting, meaning that all series of a study are retrieved completely from the archive without interrupting the load operation.

☰ Download Performance

Choose between optimal and maximum download speed.

2.7.2 Separators for Windowing Default Settings

☰ Presets Separator

Defines the separator for the window default settings.

☰ Presets Subseparator

Defines the separator within a window's default setting.

2.7.3 Modality-Depending Defaults

The individual defaults are separated by the character entered for **Presets Separator**. Within a default, the values for describing window center and width are separated by the character entered for **Presets Subseparator**.

☰ CT

Window defaults for CT series. For example: Skull,25,95].

☰ MR

Window defaults for MR series.

☰ NM

Window defaults for NM series.

PT

Window defaults for PT series.

2.7.3.1 Additional information

Window Modality Presets

The correction of the brightness and contrast of an image is specified by leveling in the DICOM standard. The window is defined by two parameters: center and width.

For most DICOM images, a window is correctly included and evaluated during image display. However, for some images, these window values are not stored at all or are poorly adjusted for the image. To correct these types of images, *syngo.share* view provides window defaults for certain modalities (called Window Modality Presets).

Format

These presets are entered in a hidden **Advanced** **Modality-depending Defaults** configuration entry and consist of a DICOM tag name followed by a regular expression for the value of the tag, the Window-Center and the WindowWidth, separated by the pipe character "|". These sequences can be repeated as needed to specify additional combinations of a DICOM tag/value, separated by a comma. If the specified tag and its contents are matched in a DICOM image, the specified window is used instead of the one contained in the image.

```
WindowModalityPresets := <WindowModalityPreset 1> ... ,  
                        <WindowModalityPreset n>  
  
WindowModalityPreset := <DICOM Tagname>|<DICOM Tag  
                        RegEx>|<WindowCenter>|<WindowWidth>
```

Example 1: Setting Window Modality Presets

```
InstitutionalDepartmentName|Trauma Room[[:digit:]]+[2911|1699],InstitutionName|General Hos-  
pital Edoras|2666|1935
```

In this example, the WC: 2911, WW: 1699 window is used for the Trauma Room department + any number. The WC: 2666, WW: 1935 window is used for the General Hospital Edoras institution.

2.7.4 External Applications

You can enter up to five external applications for opening certain file types.

Klick the **Add** button and enter the application name as well as the path to its executable file (click on **...** to call up the system browser).

Additionally, you can define command line parameters in the **Arguments** field, which *syngo.share* view will use to send certain information to the external application in question.

Hover the mouse pointer over the **Arguments** field to see a tooltip containing the most common command line parameters:

ResponseFile

Creates an RSP file containing the loaded image, including the absolute file name.


FilesDirectory

Passes the directory containing the loaded document.

FileName

Passes the file name of the loaded document.

Please see the manual of the third-party application for the form and order used for specifying the parameters.

Note

Parameters are not separated by spaces. If a space is necessary, it should be expressed with the escape sequence "%20". If a percent sign is needed, the escape sequence "%25" should be used.

Example 2: Starting an Application with Parameters

Given:

```
program.exe 0 "C:\Documents and Settings\USER\" x=2;y=3 100%
```

The aforementioned command should be entered in the following way:

```
program.exe 0 C:\Documents%20and%20Settings\USER\ x=2;y=3 100%25
```

Results:

```
argv[0]=-->program.exe<--  
argv[1]=-->0<--  
argv[2]=-->C:\Documents and Settings\USER\ <--  
argv[3]=-->x=2;y=3<--  
argv[4]=-->100%<--
```

Please read the safety notices, see [Section 1.3](#).

2.7.5 Export Application

ⓘ Path of export application

This field is used to enter the path to the export application. You can directly pass the images loaded in the print layout to an external program following export. Click on to select the path to the executable file in the file directory.

Please read the safety notices, see [Section 1.3](#).

ⓘ Show export dialog

If this option is not checked, you can directly export the print layout without opening the dialog box.

Further information about the configuration of selected export applications:

- Windows Explorer Command line parameters:

```
/e, <Responsefile> ==> opens the directory containing the document and  
creates an RSP file  
/e, <FilesDirectory> ==> opens the directory containing the document
```

Path to external application:

```
c:\windows\explorer.exe
```

- MediCAD Command line parameters:

```
/dcm <FileName>
```

Path to external application, e.g.:

```
c:\programs\medicad II\medi.exe
```

- Four Voice Command line parameters:

```
<UserLogin>&<PatientID>
```

Path to external application, e.g.:

```
c:\programs\fourvoice\fourvoice.exe
```

2.7.6 icoscript

The icoscript scripting language automates *syngo.share* view functions. It runs in the background and is used to hang images automatically and to support workflows in different departments of the hospital. icoscript cannot be disabled.

☰ Keep Layout on Hanging

Enabling this option keeps the layout when hanging images.

☰ icoscript ID

This field is used to enter the ID of the department-specific icoscript. The default script controls the automatic hanging of images.

2.8 Print

Please read the safety notices, see [Section 1.3](#).

2.8.1 System Printer

Here you can make adjustments for printing images on paper.

☰ Upper, lower, left, right margin

These values specify the gap between the page margin and the image. The default is 25 mm on each side, which is useful for A4-size pages.

☰ Resolution

This value specifies the print resolution. The higher the value, the better the image quality. Available options are 150, 300, 600, and 1200 dpi (dots per inch).

☰ Default background color

This entry allows the standard background color of printouts to be changed.

☰ Default text color

This entry allows the standard text color of printouts to be changed.

2.8.2 DICOM Printer

Set up and configure new DICOM printers and delete unneeded or unavailable printers. Follow these steps to set up a new DICOM printer:

1. Click to open the **New DICOM printer** window, enter a name and confirm with .
2. Configure the settings for the printer:

⌵ Sender AE

Application entity of *syngo.share* view. This setting is not required for all printers. Please consult the DICOM Conformance Statement of the respective printer model to determine if your printer requires this setting and how it must be configured.

⌵ Printer AE

Application entity of the respective DICOM printer. The AE is configured directly at the printer. For more information, please consult the DICOM Conformance Statement of the respective printer model.

⌵ Host

Host at the DICOM printer (may be a gateway for multiple printers).

⌵ Port

Port at the DICOM printer.

⌵ Priority

This setting is required depending on the printer. High is the default setting because certain printer models delay output otherwise.

⌵ Single-page processing

This option is designed for printing large volumes of data, e.g. slice image series with many images. If the option is enabled, the system sends only the information for one page (transparency) to the printer at a time. When processing and printing of the page are completed, the system sends the information for the next page. Single page processing prevents memory problems in case of large volumes of data.

⌵ Resolution

Please see the DICOM Conformance Statement of the respective printer model for the dpi resolution supported by the printer and enter the value. If you enter a higher value than supported by your printer, the printout takes more time because more memory is needed. If you enter a lower value than the maximum supported by your printer, the printout quality is lower than possible. The default setting is 300 dpi.

⌵ Sheet limit

If a print job exceeds the number of sheets set here, the user is asked to confirm the print job prior to printing. This confirmation is designed to prevent unintentional misprints with too many sheets.

2.8.3 DICOM Printer Settings

⌵ Interpolation algorithm

Select the interpolation algorithm used to calculate the images to be printed. The following options are available:

↳Bresenham

The default is Bresenham – it is the fastest available interpolation algorithm.

↳Bilinear

Benefit: optimum image display, but slower than Bresenham.

↳Smoothscale

Benefit: optimum image display when the scale factor of an image is below 100%.

↳As Viewer

The as Viewer setting uses the interpolation algorithm configured under  that is used for image display in the Viewer.

☰ Default background color


This entry allows you to change the standard background color for printouts.

☰ Default text color

This entry allows you to change the standard text color for printouts.

☰ Spool Directory

This is the directory where *syngo.share view* stores the individual pages of a print job if you have enabled the **Spool print jobs** option. After a page is printed, the temporary file is deleted.

Note

Caution! If *syngo.share view* crashes while printing, the temporary files are retained and must be deleted manually. The standard Microsoft Windows temp directory is selected by default. Please note that multiple gigabytes are available for spooling to prevent running out of hard drive space when printing large volumes of data.

☰ Spool Print Jobs

This lets you decide if the hard drive is used for spooling the individual pages of a print job. In this case, only additional constant main memory is needed regardless of the number of pages. If you enable this option, you must specify a directory on the hard drive as the spool directory, where the pages are temporarily stored.


2.8.4 Print Layout

☰ Allow configurable layout presets

Configure the layout buttons in the print layout. See also [Section 3.4.4](#).

☰ Automatically set focus to print layout when using camera tools

Selecting this checkbox will automatically move the print layout tab to the foreground as soon as you photograph images in the Viewer.

Note

The **Print layout** area is not moved to the foreground if you take Viewer photographs on that screen which also displays the print layout. This applies both to single-monitor workstations and the control monitor at multi-monitor workstations, e.g., radiological reporting workstation with three screens.

2.9 DICOM Export Profiles

syngo.share view supports setting or randomizing DICOM elements when exporting DICOM documents. This functionality provides what the DICOM Standard refers to as “de-identification”. DICOM export profiles provide an expedient way to change DICOM elements in repetitive exports. Basic DICOM export profiles can only be changed by administrators and can be used as templates for additional profiles. At least one base DICOM export profile exists. The export functionality itself is explained in detail in [Section 3.4.5](#).

To create a new DICOM export profile:

1. Click OR inherit rules from an existing profile by selecting one from the list and clicking .
2. Provide a unique name for the profile.
3. Select a default rule:
 - ⋮ Keep all elements
Standard Data Elements and private elements will remain intact.
 - ⋮ Remove private elements
Private Data Elements⁴ will be removed.
4. Simply type the DICOM element name or DICOM tag number to add a rule.

The following alterations are possible:

⋮ Clear

Clear field value.

⋮ Remove

Remove DICOM element.

⋮ Set random value

Randomize field value. Active only for supported DICOM elements.


⋮ Set value

Set field value to specified value. Active only for supported DICOM elements.

Administrators can define both basic DICOM export profiles and custom DICOM export profiles in *syngo.share* Configuration.

⁴Private Data Elements have an odd group number larger than 0008. These elements vary between manufacturers.

3 Usage Instructions

Note

These usage instructions provide a comprehensive description of the functions of *syngo.share* view. Due to the configuration settings of your organization or missing permissions, some of these functions may not be available to you. If you have any questions, please contact your system administrator.

3.1 Getting Started

3.1.1 User Interface

The user interface is divided into the following areas:

- [Main Menu](#)
- [Main Area](#)
- [Document Tree](#)

3.1.1.1 Main Menu

The main menu is used for navigation between the different areas and it contains the following sections, depending on permissions:

☰ Home

For logging in and out, contains release and application information

☰ Load

Main working area for loading data

☰ View

Main area for viewing and editing images

☰ Hide Patient Data

You can quickly obscure patient data in the user interface.

☰ Print

Print layout for preparing images for printing

☰ Admin

Administration area for patient merging

☰ Job Queue

Shows a list of executed actions and their status. You can clear this list by clicking on the corresponding button. The symbol in the main menu also provides a visual feedback so that you do not have to switch to the section in order to see the status of an action.

☰ Settings

This is the settings area of the application.

☰ Help

This area contains this help documentation.

☰ User

Options to log in/out, user information, application information.

3.1.1.2 Main Area

In the main area, different areas are available depending on the selected element in the main menu – e.g. the **Load** area contains various search fields and a thumbnail view for the selected results.

3.1.1.3 Document Tree

The Document Tree is visible in the areas **Load**, **View** and **Print**, and displays the loaded data. For detailed information on the Document Tree see [Section 3.2.4](#).

3.1.2 Display Modes

The following functions are permanently available in *syngo.share* view:

☰ Fullscreen

In fullscreen mode, all other components will be hidden in order to maximize the viewing area.

☰ Multiscreen

The button makes it possible to toggle between single- and multiple-screen operating modes. Multiple-screen mode automatically detects whether multiple screens are connected.

☰ White screen

This function turns all screens white, with exception of the main screen showing the view main window. This function is useful for viewing printouts of x-ray images.

3.1.3 Login and Logout

You will be automatically asked to log in, or you can manually log in via the **Home** or the **User** section. To log in, enter your user name and password, select the domain (if applicable) and click on .

If your system administrator has activated the automatic login with your Microsoft Windows user credentials, you will be automatically logged in to *syngo.share* view.

Login is not required to view locally stored data such as images from the hard drive. However, you need to be logged in to *syngo.share* view in order to view data from the archive. Always terminate logged-in sessions by clicking on to prevent unauthorized access.

At/after login, the following events may occur:

☰ Login message


You may be shown a login message. Click on to finish logging in to *syngo.share* view. If you accept the message once, you will not receive it again until your system administrator sets a new login message. If you click on , your login will be canceled.

☰ Incorrect password entry

If you enter your password incorrectly too often, you will receive a message that your user account has been locked for security reasons. Only your system administrator can unlock your account.

☰ Screen lock

After a period of inactivity, *syngo.share* view may automatically lock the screen in order to prevent unauthorized access. You can also manually lock *syngo.share* view by clicking on **Lock** in the **Home** or **User** section. To continue working, enter your password on the lock screen. This unlocks *syngo.share* view and you will see the current status of your work. You can also log out the current user on the lock screen in order to allow another registered user to log in.

Note

Note that logging out the current user means that this user's unsaved work will be lost.


3.1.4 Changing the Password

To change your password, go to the **Home** section and click on the corresponding button.

Enter the following information:

- Old password
- New password
- Confirm new password

The **Change** button will remain inactive until you have filled all fields according to the requirements.

Note

Pay attention to the password requirements. You may be required to choose a password of a certain length and use (or refrain from using) certain characters. Only if your chosen password satisfies these requirements can it be successfully changed. Passwords are case-sensitive.

Passwords can also expire if this has been configured by your system administrator. In this case, you will have to choose a new password once your old one becomes invalid.

Because passwords provide access to patient data, please remember to keep them confidential and change them regularly.

3.2 Searching and Loading

The **Load** area contains various types of searches in separate tabs. The visibility of these tabs depends on the settings of *syngo.share* view as well as your login status. The input methods as well as the presentation of the search area and the results list will vary depending on the type of search.

☰ Patient search

Use the **Archive** tab to search for documents in a patient context. This search is unavailable if you are logged out.

☰ Entity search

The entity search is analogous to the patient search. Available entities (e.g. experiments, clinical studies) are shown by name in separate tabs, provided that the entity search feature has been enabled in the settings. The appearance of the search area varies depending on the entity in question and comes with its own predefined search fields. However, you can still adjust these search fields manually.

☰ Collections

In this tab, you can search for document collections created for specific purposes, e.g. morning briefings, using a variety of criteria such as the creator, their modification date, or share settings. For more information on collections see [Section 3.2.4.5](#).

☰ DICOM Q/R

In this tab, you can search an external system (external PACS), provided that one has been configured. A DICOM Q/R search is conducted similar to a patient search, however, its results list is completely independent of the patient search results list.

☰ External search

In the **External** tab, you can start a cross-enterprise search (XDS search) for external patients and load external documents. This search is only available if the eHealth Solutions Consumer Adaptor has been configured in *syngo.share* view.

☰ Filesystem

In this tab, you can search for files available on your local workstation (on hard drives, flash drives, etc.). You can also search for DICOM files (e.g. on a patient CD), see [Section 3.2.2](#).

3.2.1 Searching in Archives

The patient, entity, collections, DICOM Q/ R and external searches share certain features and input methods.

In order to start a search, enter the desired criteria (the input varies depending on the search field type, see [Section 3.2.1.2](#)). Use next to a search field in order to reset it or press in order to reset all search fields.


Click on or press to start your search.

Keep in mind that the entry of multiple search criteria will naturally lead to a narrower search result. You can also save your search query, see [Section 3.2.1.4](#).

The results will be shown in list form beneath the search field area.

The results area is composed of **Studies** and **Series**. The **Studies** area contains the actual search results list. Found entries in the results list can be DICOM studies (examinations) or generic documents (generic containers). For more information, see [Section 3.2.1.5](#).

The **Series** area offers a preview and detailed information about the elements within a selected DICOM study or container. You can adjust the area size by clicking and dragging the **☰** symbol. For more information on the **Series** area, see [Section 3.2.1.6](#).

Note

If you cannot find patient data which you know exists and you have entered the search criteria correctly, the examination you are looking for may have been performed outside the organizational unit you belong to. In this case, please ask your system administrator for authorization to view data from other organizational units.

3.2.1.1 Configuring the Search Fields

syngo.share view offers a multitude of search fields which you can freely customize according to your needs, e.g. by removing those fields which you rarely use or by changing the order of appearance to prioritize your preferred search fields.

To customize your search fields, click the arrow symbol next to **Search** and select **Configure search fields**. This will open a dialog box containing all available and all displayed search fields. Use the arrow buttons to add/remove displayed fields and rearrange their order of appearance.

Note

Your system administrator may have designated certain search fields as mandatory. These fields will be highlighted.

3.2.1.2 Types of Search Fields


Depending on the category a search field belongs to, input methods vary and certain additional conditions may apply.

Free Text

You can enter search terms directly into this type of field (e.g. First and Last Name, Referring Physician, Patient ID). The search is case-insensitive. You can also use wildcard characters to widen or narrow the scope of your search.

↳ Wildcards % and *

The % or * character can be used as a placeholder for any number of characters and is always automatically (and invisibly) appended to the search. For example, entering the letter “m” in a free text search field will list all entries starting with this letter. However, a search using only % or * is not possible.

Note

This automatic addition of % does not occur in ID fields such as Patient ID, Accession Number or Visit ID in order to avoid false positives when searching for patients via barcode scanner.

Example 3: Search using %

If you search for last names using “W%r”, the results will include not only names like Wheeler and Walker, but also Waterfield, Whorton, Workowski, etc.

↳ Wildcards _ and ?

_ or ? act as a placeholder for a single character, which allows for a narrower search than % and *.

Example 4: Search using _

If you search for last names using “m_ller”, you will receive results for e.g. Miller and Muller, but not Mueller. However, longer names such as Mellerson and Mullerbach are still included due to % being automatically appended to the search.

Special Free Text Search Fields

↳ Mother’s Maiden Name

The mother’s maiden name may contain special characters. Due to technical reasons, the search field does not support some of these special characters. If the name you are searching for contains one or more of the special characters `\ | ~ & ^`, please use a wildcard at this position to conduct your search.

↳Description

This field searches the DICOM tag StudyDescription. The % wildcard is automatically prepended and appended to the search term.

↳Patient Age/Study Age

These fields require a particular syntax. In addition to entering the desired number, you have to indicate the desired unit of time (y or a = years, w = weeks, d = days, h = hours). Searches with only the numerical values are not possible. By adding > or <, you can search for certain age ranges. Without a > or <, the search will only include the precise age.

Example 5:

Searching for patients younger than 20 weeks: <20w

Searching for patients exactly 70 years of age: 70a or 70y

☰ Selection List

Fields of this type (e.g. Sex or Document Type) require you to choose an entry from a predefined list. Certain selection lists contain items with subcategories. Simply click the arrow symbol next to the item in question in order to view and select a subcategory.

Certain selection fields, e.g. Modality or Organizational Unit, allow multiple list entries to be selected in order to expand the search.


Special Selection List Fields

↳Region

This search field filters for the body part or body region examined. The available items are based on the DICOM tag BodyPartExamined.

☰ Date

Date fields offer several options: You can either set the current date by simply clicking in the field, input the date manually, or select a date from the calendar by clicking on the arrow symbol.

Note

The date format and calendar entries depend on the format and region configurations of the operating system on the local machine (the Microsoft Windows configuration option is found under **Control Panel** > **Clock Language and Region** > **Date and Time** > **Region**).

Special Date Fields

↳Study Date

This field combines a selection list with a date search. First, select a time period indicator from the list. If the indicator is connected to the current date or time (e.g. "today" or "yesterday"), you can start the search immediately. However, if you wish to select a certain date ("on", "from", "until") or a certain date range ("between"), one or two date fields will open where you can enter the desired date as described above.

3.2.1.3 Searching by Keyword

You can also search for documents according to their keyword tags. The configurable keyword fields depend on the available keyword groups.

To search for a keyword click on the arrow button in the keyword search field. A list of keywords from all associated keyword catalogs opens. To filter the list, enter part of the desired keyword. Click on the desired keyword to add it to the search field (you can add more than one keyword to a search). Keyword search fields associated to free text keywords do not provide a selection list, you have to enter the search terms manually.

A search for keywords from catalogs is always an exact search. For searches in a free text catalog however, a placeholder is always added automatically at the beginning and the end of the text entered into the search field. These search terms are set in italics in the search field.

All documents will be found that are tagged with at least one of the specified search terms on at least one tagging level, i.e. if e.g. one instance of a series is tagged with a keyword searched for, the entire corresponding study will be included in the results.

To remove a particular keyword, click the right next to it.

For more information on keywords, see [Section 3.2.4.4](#).


3.2.1.4 Saving a Search Query

You can save frequently used search queries to avoid having to reenter them every time. Doing so will also save the entered search terms, which will be restored when the saved search query is selected.

Example 6: Display all CTs of a given day in the results list

If you have a daily CT meeting and want to see all CTs of that particular day in the results list, you can create a CT_Meeting search query, save it, and retrieve it by clicking on its name.

1. Fill in the desired fields in the search area.
2. Click to start the query.
3. If you are satisfied with the result, you can click on to save it.
4. Enter a name for the query (max. 50 characters, comma and semicolon are invalid) and click .
5. You can find the saved search query as a new button above the search field area. Clicking on it will execute the saved query.

Note


Search queries in entities, collections and DICOM Q/R are saved locally on your workstation and are available to all users of the same computer.

Setting a Search Query as Standard

You can set a saved search query to start automatically as soon as you log in to *syngo.share* view. Click on the arrow next to the desired search query button and select from the submenu. The standard search will be starred. At your next login, *syngo.share* view automatically switches to the appropriate tab and loads the results list of the standard search query.

Editing a Saved Search Query

Click on the arrow next to the button of a saved search query to modify it. You can edit the displayed fields, rename the query or delete it.

**Note**


As long as you are logged in, *syngo.share* view saves your searches and results lists automatically. You can switch between them using the arrow buttons on the left above the search field area. When you log out of *syngo.share* view, these temporary search queries are deleted.

3.2.1.5 Results List


Similar to the search fields, you can configure the results list according to your needs. By clicking on , you can add or remove columns from the results list or change their order of appearance. For detailed information on configuring the displayed items, see [Section 3.2.1.1](#).

In addition, you can

- drag columns to change the order,
- adjust the column width manually by dragging their border or automatically by double-clicking on it,
- sort the columns by clicking on them (the displayed up- or down-arrow indicates the type of sorting in use).

**Note**

Keep in mind that the results list is generally static, meaning that it does not reflect changes to archive data that are made after a search query has been sent. If you wish to update your results list, either repeat the search manually or set an update interval (see [Section 3.2.1.5](#)).

**Note**

Due to visibility restrictions you might not be authorized to view all parts of a document. This might lead to an adjustment in the series and/or object count once you click on an entry in the results list.

Search Results List Context Menu

The context menu offers a means to further organize and filter the search results list. It can vary according to your permissions, the type of search and the selected document type.

Grouping

You can group the results based on up to five columns. These columns have to be arranged in the exact order according to which grouping is to be performed. Grouping will subsume entries with matching values within one column into a group. In case of multiple groups, grouping is performed hierarchically from left to right. The grouping icon  appears in the header of every grouped column, indicating the number of grouping levels.

Search for prior examinations (Patient)

Searches for all prior examinations of the patient currently selected in the results list.


Search for examinations (Patient)

Searches for all examinations of the selected patient.

Search for prior examinations (Visit)

The Visit ID is used to search for prior examinations belonging to the same case as the current examination.

☰ QR Code

Displays the patient QR code. This option is only available in the **Archive Search** and **DICOM Q/R** tabs. You can scan the QR code with the mobile app Visual Capture to quickly find the correct patient and thus add health data gathered with your mobile device (image, video and audio files) to their health record. You can also display the QR code by clicking on the  symbol next to a patient entry in the Document Tree or by selecting the option in the context menu (see also [Section 3.2.4](#)).

☰ Workflow State

Here you can change the workflow state of a document, provided you have the appropriate permission. The available workflow steps are defined by the system administrator in *syngo.share* Configuration.

Example 7: Possible Steps of a Workflow

- New
- Prepared
- Reported
- Signed-Off


☰ Tags

You can add keyword tags to documents. For further information, see [Section 3.2.4.4](#).

Copying Results to Microsoft Excel

You can export search results to Microsoft Excel by selecting the desired items from the results list. Press **Shift** while clicking lets you select multiple consecutive entries; pressing **Ctrl** while clicking allows you to select multiple non-consecutive entries.

If you are satisfied with your selection, either use **Ctrl + C** or click on  to copy the items to the clipboard.

Note

Please note that the context menu entry **Copy** will only copy the contents of the results field from which the menu was called.

Now the copied items are ready to be inserted into a new Excel worksheet by clicking on a table cell and using **Ctrl + V** or by clicking the paste icon in the toolbar.

Setting Automatic Updates for the Results List


You can make the static results list more dynamic by setting an automatic update interval. Click on  and choose one of the available intervals.

The set interval applies to saved search queries as well as manual searches. If a saved search query has been set as standard, the update interval starts when you log in to *syngo.share* view.

Warning

The automatic update of the results list causes a higher system load. The shorter the update interval, the higher the load. To avoid an unnecessarily high system load, choose longer intervals or, if unneeded, disable the automatic update function altogether.

Note

If the automatic update is deactivated, changes to the data are not shown between searches. If you wish to update the results list manually, simply repeat the search by clicking .

Restricted Document Parts

Documents may contain restricted elements, i.e. series or generic files which you do not have permission to view. These restricted elements will not be shown in the preview area or the document tree, and cannot be loaded for viewing. However, depending on your other permissions, you can still perform operations which affect the entire document, such as copying, moving or sharing the document.

Documents with restricted elements are indicated in the **Objects** column of the results list, which lists the total number of elements contained in a document (series or generic files). When you select a document with restricted elements, its total number of elements will be divided into visible and hidden elements.

3.2.1.6 Using the Series Area in the Archive/Entity/Collections Search

The series area is available in several search tabs to display details for selected results. Depending on the type of search, its layout and the information available may vary.

In the archive, entity and collections searches, the series area offers a thumbnail view of the files contained within the selected DICOM studies, generic containers or collections. These thumbnails are generated when the data is added to the archive. Preview images of files (provided they are available) can only be displayed for files that can be viewed and edited in *syngo.share* view. Files from external applications are displayed using their generic file icon.

File Information

The thumbnail view offers a lot of useful information on the selected files before loading them into the Document Tree. (Files in the thumbnail view do not count as loaded.)

≡ Image number/Filename/File count

Depending on the presentation options, files in the thumbnail view are displayed with their image number (DICOM series) or name (generic document) or the number of files in a stack (only available in stack view).

≡ Stackability

The stack symbol indicates elements that are stackable or, if stack view is activated, elements that have been stacked. Note that generic files cannot be stacked.

≡ Tagging

An orange triangle lets you know which files have been tagged with keywords.

≡ Information on hover

By pointing your mouse at the desired element, you can display a tooltip containing relevant information such as patient ID, study date, document description, etc. The information available depends on the selected element.


☰ Version

Hovering the mouse pointer over an element will also display an ☰ symbol. Click on it to receive information about the document version.

☰ Quickview


Hover the mouse pointer over an element to display an **Q** symbol. By clicking on it, you can open a quickview window which offers a slightly larger preview image, detailed information on the element as well as rudimentary viewing tools. If stack view is activated, you can scroll between elements in a stack by pointing at the slider bar and scrolling with your mouse wheel.

If you open several images one after another in quickview, the most recent image is displayed on top, hiding all the others. Simply move the quickview window(s) in order to view the images side by side. However, a more convenient method for side-by-side viewing is to open the images in the Viewer and activating the comparison mode. See [Section 3.3.1.4](#) for details.

Note

Quickview offers the same viewing and editing tools that are also available in the Viewer. Calculations done in quickview are not saved automatically; however, they are retained if you load the associated documents in the document tree.

Note

If quickview cannot display a particular generic file (e.g. a Word file), clicking on **Q** will open the file in an appropriate third-party program, provided it is installed on your PC.

Display Options

You can customize the series area according to your needs.

☰ Hide series area


You can hide the series area using **F8** or by clicking on the ☰ symbol.

☰ Filter according to ICOTAGS keywords

You can filter the elements shown in thumbnail view by choosing to display only those tagged with a certain keyword from the ICOTAGS catalog, e.g. "findings-relevant". Simply choose the desired keyword from the selection list. For more information on keywords and ICOTAGS, see [Section 3.2.4.4](#).

☰ Grouping and sorting

Grouping documents is useful when e.g. multiple series belonging to different patients are shown in the thumbnail view. Select the desired entry from the grouping list to group the thumbnail elements beneath their respective headers (e.g. grouping elements by patient arranges the elements beneath headers bearing their associated patient's name). You can also change the sorting of a group by ticking the checkbox **Sort descending** (the checkbox is deactivated by default, meaning that groups are sorted in ascending order).

Note

Grouping according to study always shows the most current study first. The older the study, the further down it is found on the list. The same is true for grouping according to examination date. If the **Sort descending** checkbox is activated, older groups are shown first.

Stack series

You can use the **Stack series** checkbox to achieve a better overview of series containing a great number of images (>100). Stackable elements are shown with the stack symbol in the thumbnail view. When stack view is activated, the number of files contained in a stack is displayed (e.g. #167) alongside the stack symbol.

Adjust thumbnail size

To increase or decrease the size of the thumbnails, either drag the slider to the left or right or click on the  and  icons.

3.2.1.7 Using the Series Area in the DICOM Q/R and External Search

The series area in these two search tabs differs from the regular thumbnail view. If you select a study in the DICOM Q/R search results list, the series area shows detailed information on the series contained therein in several columns. You can adjust the displayed columns just like the columns of the search results list (see [Section 3.2.1.5](#)). Similar to the thumbnail view, elements shown in the DICOM Q/R series area are not loaded.

3.2.2 Searching Local Storage Media

The local filesystem search differs substantially from the archive searches. Instead of configurable search fields, the filesystem search offers a file browser which functions just like Windows Explorer.

The left-hand area contains the file browser for navigating the visible drives and folders while the right-hand area shows the elements contained in the selected device or folder.


You can select a directory to browse in one of two ways:

- Scroll through the navigation tree and click on the desired folder.
- Click on  in the header bar to quickly select a specific folder, device or the Windows Favorites folder.

You can set a selected directory as your start directory by clicking on .

Use  to update the filesystem browser display in order to see any changes to the folder contents that may be happening in the background, e.g. the deletion of a file.

The area below the file browser offers means of adjusting the layout of the right-hand area. You can switch between list view and icon view. The list view offers detailed information on the individual elements via sortable columns. You can drag the slider bar or click on  and  to adjust the size of the preview images. The DICOM overview (to be selected via the CD symbol) shows all DICOM files in the selected directory.

Just like in Windows Explorer, you can select multiple files with the help of  or . The area below the file browser also offers a  button.

Just like other types of searches, you have to load the desired files into the Document Tree in order to view and edit them (see [Section 3.2.3](#)).

Searching for DICOM Files (e.g. on a Patient CD)

You can also search a selected directory or medium (e.g. a patient CD) for DICOM files by clicking on the CD symbol in the footer bar (). If the selected directory or medium contains a DICOMDIR directory structure, all DICOM files included therein are listed immediately. In case no DICOMDIR directory file is found, you can scan a CD, USB stick or any directory manually.


3.2.3 Loading Search Results

Before documents can be opened and edited, they have to be loaded in the Document Tree and hung. You can quick-preview documents in the series area or file browser without loading them in the Document Tree first.

Depending on the type of search, you have various options for loading results into the Document Tree.


When searching local storage media, you can load elements by double-clicking, via drag-and-drop, or by using the available buttons below the browser area.

When searching in an archive (patient search, entity search, DICOM Q/R, etc.), you can load elements directly from the results list or from the series/thumbnail view area. To load a desired element, either double-click it, use the shortcut **F9** or drag the result to the Document Tree using the mouse. Using the latter method, it is possible to select and load several results at once. All loaded documents are shown in the Document Tree ([Section 3.2.4](#)).

Note

If not all images of a series are loaded, a threshold for the number of images to load was set in **Settings > Advanced**; see [Section 2.7.1](#).

Note

If the result is listed as offline (indicated by the yellow dot in the results list status column as well as in the thumbnail view), you are asked if you want to initiate the retrieval. If you confirm, the documents are retrieved from the long-term archive, which could take some time, as will be indicated in the Document Tree. As soon as the examination is fully available, it will be displayed with its corresponding thumbnail preview and can be hung and edited as usual. You can check the retrieval status in the **Job Queue** area in the main menu.

3.2.4 Document Tree

The Document Tree on the far right side displays all loaded elements. DICOM documents are highlighted in blue whereas generic documents are highlighted in green.


The display hierarchy for DICOM documents is always the same:

Patient > Examination > Series or Document

For generic documents, this hierarchy may differ.

The individual elements are displayed with a preview image (if available) or a generic file icon alongside additional information (e.g. patient ID, birthdate). The additional information shown depends on the settings made in **Settings > Labels > Loaded Documents** and can be adjusted if necessary.

A progress bar shows the loading progress. If necessary, you can halt the loading operation by clicking on **Stop**. A stopped loading operation can be resumed by clicking on **Play**.

Note

If a download threshold has been set in **Settings > Advanced > Memory Management**, only the configured number of images per series will be loaded. Click **▶** to manually load additional files.

The **🔑** symbol indicates that a document was tagged with keywords. Clicking on the symbol will open a window where you can view the assigned keywords. For more information on using keywords, see [Section 3.2.4.4](#).

3.2.4.1 Basic Actions

You can perform the following basic actions in the Document Tree:

☰ Open element

To open and view an item from the Document Tree, simply double-click on it. For special functions, please consult the following section. To learn more about hanging documents, refer to [Section 3.3.1](#).

☰ Close element

To close a single image or document, a series, a study or a patient entry and remove it from the Document Tree, click the  on the appropriate level.

☰ Empty Document Tree

In order to close all loaded elements in the Document Tree, simply click the big  in the Document Tree menu or use .

3.2.4.2 Document Tree Main Menu

The main menu is displayed above the Document Tree itself and offers various presets for loading and viewing documents in the Viewer, as well as functions for elements in the Document Tree:

☰ Layouts

These buttons can be used for quickly defining the layouts available in the Viewer. You can choose from various pre-defined layouts or define one yourself.

☰ Action before loading

This selection list determines if previously loaded documents will be removed before loading new documents into the Document Tree.

☰ Action after loading

You can determine whether a certain action should be performed after documents have been loaded into the document tree. For more information see [Section 3.3.1](#).

☰ Patient selection

If several patients are listed in the Document Tree, you can filter them using this selection box. This means that only the documents associated with the selected patient are shown in the Document Tree. Using the arrow keys, you can quickly navigate the list.

☰ Sorting

You can sort loaded documents via  **Configure**. Sorting is possible at the studies, series and files level. The main sorting criteria are field type and document properties, which can be further defined using the second selection box. If the **Descending** checkbox is ticked, the items are sorted in descending order.

☰ Save collection

 **Save collection** saves the active tab in the Document Tree as a collection and/ or allows you to save changes made to an existing collection. For more information on collections, see [Section 3.2.4.5](#).

☰ Export functions

 offers several export functions, e.g. the export of selected series. For more information, see [Section 3.4.5](#).

☰ QR code

You can show the patient QR code on the patient level by clicking on  in the Document Tree main menu or selecting it in the context menu.

3.2.4.3 Document Tree Context Menu


The available entries in the context menu may differ depending on the selected level and document type.

☰ Search for Examinations (Patient)

On the patient level, you can directly search for all other examinations associated with a particular patient and subsequently load these results.


☰ Search for Prior Examinations

On the patient level, you have several options for searching and loading suitable prior examinations to a current examination (study). The variant entries of [Search for prior examinations](#) can assist you in refining your search.

Note

If you wish to compare prior examinations to the current examination, use [Compare progression](#) or [Compare with current](#) on the series level. For more information, see below.

Note

The search for prior examinations does not use the **Description** criterion. For example, if the description for the current examination is KNEE, the system does not only search for this particular patient's prior knee examinations. Since descriptions can vary widely, such as KNEE or KN, THORAX or THX, the system may be unable to find prior examinations with different descriptions. If necessary, you can specify your search by adding additional criteria, such as Document type = DICOM instead of dcm to search only for DICOM studies, or by using a precise description such as KNEE.

☰ Hang

This command is available on all levels of the Document Tree and is used to open either a document, series, study or all loaded documents associated with a patient in the **Viewer**. For more information see [Section 3.3.1](#).

☰ Series overview

This command is available on all levels and opens series in the Viewer using overview mode (see [Section 3.3.1](#)).

☰ Change description

The document description and study date can be changed for studies or generic documents.

☰ Delete

A study, individual series of a study or generic documents can be deleted from the archive. Since this deletion is permanent, you will receive a prompt to confirm your choice.

☰ Send

A study, series or generic documents can be sent directly from the Document Tree (see [Section 3.4.9](#)).

☰ Tag

Tagging is possible on study or series level, or for generic documents. The keywords are directly applied and saved. For more information on tagging, see [Section 3.2.4.4](#).

☰ Compare progression

Show documents on multiple monitors in their series order. For more information on comparing images, see [Section 3.3.1.4](#).

☰ Compare with current

Keep the current document on the first monitor. The next monitor shows the series of an earlier examination for comparison.

☰ Compare mode

This mode will show the series of the latest examination on the first monitor. The next monitor can display the series in different ways.

3.2.4.4 Keyword Tagging

syngo.share view makes it possible to apply keywords to documents. This allows you to search for and group documents by their keywords.

Keywords can be based on catalogs (e.g. ICD 10) and/or on free text. The advantage of catalogs is that the keywords are precisely defined, improving searchability. The advantage of free text keywords is that you can freely add keywords as desired, though this brings with it the disadvantage of making entries more difficult to find due to the possibility of multiple synonymous keyword entries (e.g. "femoral fracture" and "broken thigh bone").

☰ Catalog


This contains predefined keywords.

☰ ICOTAGS

This special catalog allows only one single keyword to be set. This exclusive behavior is useful for things like "diagnosis-relevant" keywords.

☰ Free Text

These keywords are not bound by a predefined catalog. Any desired keyword can be used to tag a document.

Note

The ability to index documents and add keywords depends on your permissions in *syngo.share* view. If you do not have permission to do so, please contact your system administrator.


You can assign any number of keywords to studies, series and images as well as generic documents with any level from the keyword trees. Keywords can be removed as well, e.g. if they were incorrectly assigned.

To assign keywords on different levels, right-click the element on the desired level in the Document Tree and select the tag entry from the context menu. If you want to tag single images of a series, load them in the Viewer, open their context menu and select **Tag**. In both cases, the tagging input window opens, allowing you to select and assign keywords.

Keyword groups can be selected from a selection menu. If a catalog is expandable, you can enter a new keyword in the text field. It is also possible that a field has no keyword catalog assigned to it. In this case, the catalog is a free text catalog and the tagging is done by entering the keywords manually.

Using the arrow button next to the fields, the respective keyword catalog of a keyword group is opened. Filtering allows you to search for keywords in the catalog. Tagging can be done by double-clicking on an entry in the catalog or by clicking **Apply**. Clicking **OK** saves your tagging.

To remove a keyword, click the info button next to an element in the Document Tree, or click the orange triangle in the Viewer. The tagging window appears and keywords can be removed by clicking **x** next to the keyword.

Note

The assignment or removal of keywords occurs immediately, meaning you do not have to save it explicitly.

For information on keyword searches, see [Section 3.2.1.3](#).


3.2.4.5 Collections

These are collections of documents that can be prepared and saved, e.g. for morning briefings, and can be retrieved and opened at a later date.

☰ Create Collection

Collections are prepared in the Document Tree:

1. Click on **+** in the Document Tree, enter a name for your new collection and confirm with **OK**.
2. Your new collection now appears as a separate tab in the Document Tree.
3. Next, search for all examinations you want to add to the collection and open them in the Document Tree.
4. Save the new collection by clicking on **☰** and choosing **Save collection**.
5. A new window opens where you can edit the name and add a description (max. 800 characters).
6. You can choose to make your collection available to other users and/ or user groups. Click on **Share with** to select individual users or click on **Share with groups** to select user groups. You can press **Ctrl** or **Shift** to select multiple users/groups at once. If you choose not to share your collection, it will remain private.

Note

Keep in mind that the newly created collection is automatically active in the Document Tree (indicated by the appearance of the tab). This means that any documents you load from the results list are automatically added to this collection. If you wish to add documents to a different collection or load results independent of a collection, you first have to switch to the desired tab.


☰ Search for and Load Collections

Use the **Collections** tab in the **Load** area. For general information on searching, see [Section 3.2](#).


Either double-click on the desired collection, choose **Load** from the context menu or use the **Open** button to load the collection in a tab in the Document Tree. Each collection is opened in a separate tab.

☰ Edit a Collection

If the collection is not shown as a tab in the Document Tree, you have to search for and load it first.

You can now execute further searches in the **Load** area and add the desired results to your collection. Save your changes via  **Change**. You can also adjust the name, description and viewing permissions.


It is also possible to keep the original collection intact and save your changes as a new collection instead. Simply activate the checkbox **Save as new collection** in dialog window.

Note

Only the owners of a collection are allowed to change the individual permissions.

☰ Delete Collection

Search for the desired collection, right-click on the result and select  from the context menu. You will receive a confirmation prompt. Click  to permanently delete the collection.

3.3 Image Display and Editing

3.3.1 Hanging Documents

Hanging images from series or documents allows you to view, compare and edit them in various ways (e.g. by using annotations). The Viewer offers predefined layouts as well as a number of display modes (band mode, compare mode, and snake mode).

To hang and view a series or document, double-click it in the Document Tree or use the context menu in this area.


The context menu offers several options for hanging documents:

- Hang
- Hang 2x2
- Series Overview
- Compare Mode
- Hang Structured Report
- Action before Loading
- Action after Loading
- Applying a Hanging Workflow

3.3.1.1 Hang


The  command in the Document Tree context menu is only available for DICOM files. For generic documents, choose .

The  command is available on all three levels of the Document Tree (patient, examination and series). Depending on where the context menu is called, the command either hangs all documents associated with a patient, the selected examination or the selected series.

**Note**

With the exception of series, **Hang** always hangs to all monitors and enables band mode.

You can also select multiple series or examinations on the respective level and use **Hang** to hang them all at once. The series are hung from top to bottom in the order listed in the Document Tree. This is the preset default method called band mode.

**Note**


Please note that a mixed selection of series, examinations, and patients cannot be hung.

Above the Document Tree, you can choose between several hanging arrangements. Symbols at the top of the Viewer area also provide options for varying the display order.

3.3.1.2 Hang 2x2

The **Hang 2x2** command in the context menu hangs a loaded examination and divides the Viewer window into 2x2 series per screen, allowing you to view four series on each screen side by side.

Just like normal hanging, **Hang 2x2** hangs all series for an examination to all monitors and enables band mode.

**Note**

The **Hang 2x2** command is shown in the context menu for loaded examinations (identified by the folder icon) only, not in the context menu for loaded series or patients.


3.3.1.3 Series Overview

This context menu item hangs all loaded series. A layout is set automatically for every series, allowing you to view as many images for the series as possible (a maximum of 4x4 images).

3.3.1.4 Compare Mode

You can activate compare mode by selecting the context menu item of the same name (this sets a checkmark).

This mode shows the series of the most recent examination on the first monitor. The series can be displayed differently on the next monitor.

**Note**

Compare mode can be disabled as needed by opening the context menu in the Document Tree and unchecking the **Compare mode** checkbox. The linking of both Viewer windows is disabled.

The following prerequisites have to be met to use compare mode:

- To compare images, you need at least two screens.

- The desired examinations/series must be present in the Document Tree.
- The examinations/series to be compared must belong to a single patient.

Compare with Current

This method assumes that only series with identical study and series description (e.g. CR Thorax Exam) are compared.

In the Document Tree, select the series you wish to compare with others and choose **Compare with current** in the context menu.

All series with identical study and series descriptions are hung from left to right chronologically by creation date, starting with the most current series.

The most current series is pinned to the first monitor. The second monitor will show the series which was used to initiate **Compare with current**.

By scrolling, you can now compare the most current series on the first viewing monitor with the matching series from prior examinations on the second monitor.

Example 8: Current Comparison

You have loaded four thorax examinations with different examination dates for a patient. Each examination contains a comparable series with a thorax image (identical series description = THORAX). **Compare with current** will hang the four series. The most current series is hung on the first viewing monitor. The therefore comparable series from earlier examinations are hung in chronological order on the second viewing monitor. When scrolling, the current series remains pinned to Monitor 1. Monitor 2 scrolls back and forth between the comparison series.

Compare Progression

This context menu command hangs all loaded series associated with a patient which have the same series and study description.

The series are scrolled through chronologically on all monitors. No series is pinned to the first monitor.

In order to use **Compare progression**, select the series you want to compare to others in the Document Tree and choose the command from the context menu.

All series with identical study and series descriptions are hung from left to right chronologically by creation date, starting with the most current series.

The first monitor will initially show the series which was used to initiate **Compare progression**.

Example 9: Compare Progression

You have loaded four thorax examinations with different examination dates for a patient. Each examination contains a comparable series with a thorax image (identical series description = THORAX). Use Compare progression to hang the series. When you scroll, the system always scrolls by one series on both screens. This allows you to view two series acquired in a chronological sequence on Monitors 1 and 2.

Compare Examinations

To compare examinations, load a current and matching prior examination into the Document Tree and select **Compare examinations** in the context menu at the examination level.

The examination used to initiate **Compare examination** is considered to be the current examination. All series from this examination are hung on the first monitor.

The matching prior examination is the prior examination which is chronologically closest to and whose series description matches that of the current examination. The series of the matching prior examination are hung on the second monitor.


Note

If you receive the message “No prior examination loaded”, you have selected an examination as the current examination for which no prior examination is available in the Document Tree. Please check if any prior examinations exist for the examination in question.

3.3.1.5 Hang Structured Report

DICOM Structured Reports are opened in the Viewer. If available, you will see links to the individual documents. You can click a link to show an image in the quickview. You can edit images as usual.

Moreover, you can view the metadata of a Structured Report or its linked documents. See [Section 3.3.7](#) for more information.

Note

If you open a Structured Report from the thumbnail view without loading it into the Document Tree, the linked documents cannot be presented.

3.3.1.6 Action before Loading

If you select **Remove documents**, every loaded patient in the Document Tree will be deleted upon loading a new one.

3.3.1.7 Action after Loading

You can choose between **Clear and hang** and **Hang patient's documents** in the drop-down menu.

3.3.1.8 Applying a Hanging Workflow

As already mentioned in [Section 2](#), you can create Hanging Workflows to execute sets of predefined actions when preset conditions apply.

Available workflows are shown in the context menu of the Document Tree as long as at least one of their preset conditions applies to the loaded document(s). For example, if a certain workflow has been set to only apply to DICOM-type OP documents, it only appears in the context menu if documents of this type have been loaded in the Document Tree.


Select the desired workflow in the context menu to automatically execute its steps.

For more details on the creation of Hanging Workflows see [Section 2.4](#).

3.3.2 The Viewer

The **View** area is the screen area used for viewing and editing (e.g. diagnosis, calculation) both DICOM and non-DICOM images.

To move the Viewer to the foreground, go to the **View** area or move items to the area by simply dragging them (only possible in multiscreen mode, see [Section 3.1.2](#)).

Note

Every image shows custom text information about the patient, physician, hospital and department, examination, admission date, and technical image data. If text information does not appear, you can enable it using the **T** key or by selecting **Viewing properties** > **Show text** from the context menu. The image labeling can also be configured individually. For more information, see [Section 2.5](#).

The Viewer can be spread across multiple monitors and also be used in fullscreen mode. To display multiple series on a screen, you can divide each screen into multiple Viewers using layouts. You can also define the layout within a Viewer with Viewer layouts, see [Section 3.3.3.9](#).

The Viewer provides different viewing options depending on the document type shown (e.g. metadata, ECG, movie) which in turn use different controls. In addition, the context menu of the Viewer is always available. For more information, please see the respective chapters.

In the **View** area, you can either place Viewers next to each other or define a layout within a Viewer. The first option is useful for viewing different series or images in different Viewers, where specialized viewing tools are available (see [Section 3.3.1](#) for a detailed description). You can use layouts within a Viewer to view the images of one series all at once, see [Section 3.3.3.9](#).

3.3.2.1 Dividing the Viewer Vertically and/or Horizontally

By clicking on the layout buttons in the Document Tree or at the top of the Viewer area, you can apply various predefined layouts to the Viewer. You can either drag and drop items from the Document Tree into the individual Viewers, or select several series or images in the Document Tree and then apply a predefined layout. The layout containing the selected documents will be opened in the Viewer.

It is possible to adjust the predefined layouts to your needs. Click on the  icon or the down arrow next to the predefined layout buttons and design your desired layout by clicking and dragging the mouse (max. 5x3).

If you have already loaded a study to the Viewer and want to change the layout, the remaining series in the Document Tree are automatically loaded and presented in the Viewer.

Using  and  above the Viewer, you can browse between the layouts. The  button clears the layouts as well as the Viewer (see also [Clear Viewer](#)).

You can move individual Viewers on a monitor or between the individual monitors. If changes such as zoom, window/level, etc. were previously applied, these changes remain intact. To move a Viewer, click its corresponding toolbar and drag the series to the desired location.

3.3.2.2 Defining the Layout of Individual Viewers

You can divide individual Viewer windows in up to 6x6 sections to view more images of a series side by side. There are several options for applying a layout:

- the context menu item **Layout**,
-  **Layouts** in the toolbox (only available if the toolbox has been configured, see [Section 3.3.2.5](#)),
- or  in the toolbar (only available if the toolbar has been customized, see [Section 3.3.2.4](#)).

Choose the desired number of fields by pressing and holding the mouse button.

To fill the Viewer with a single image, double-click the desired image to enlarge it. To return to the divided Viewer, double-click the enlarged image again.

3.3.2.3 Viewer Context Menu

Open the Viewer context menu with a right click. The appearance of the menu and the available points may vary depending on your user permissions.

The toolbox is an alternative to the context menu. Although it differs in appearance, you can access it in the same way and it contains the same functions as the context menu. For more information, see [Section 3.3.2.5](#).


Note

Many context menu entries have a key or key combination listed to the right of the text. Instead of opening the context menu and searching for a menu item, you can simply press the corresponding key (combination) in the Viewer to perform the same action as the menu item. See [Section 3.5](#).

3.3.2.4 Toolbar

The toolbar is shown to the right and below the Viewer window. To personalize the toolbar, right-click on it and choose the desired tools from the selection window by using the arrow buttons. You can also use the arrow buttons to change the order of appearance for the individual tools.

3.3.2.5 Toolbox

The toolbox is a freely configurable alternative to the Viewer context menu and is also opened by right-clicking. To use the toolbox instead of the context menu, you have to first enable it in    . These steps have to be taken regardless of whether you are logged out (at the host level) or logged in (at the user level).

The various toolbox tabs contain a practical set of tools based on the primary imaging device. Depending on the given image, the correct tab will automatically be moved to the foreground when opening the toolbox.

Some buttons have a down arrow icon which indicates that you can choose the desired tool from a sub-menu. Once you select a tool, the toolbox closes automatically.

Tabs

By default, the toolbox contains five preset tabs:

☰ Default

Automatically selected if the type of imaging device is not available in the image data.

☰ X-ray

Automatically selected for the following imaging devices:

- Computed Radiography (CR)
- Digital Radiography (DX)
- Direct Radiography (DR)
- Radio Fluoroscopy (RF)

☰ CT / MR

Automatically selected for the following imaging devices:

- Computed Tomography (CT)

- Magnetic Resonance Imaging (MR)

Special

Mainly contains menu items used to change the image (pseudo colors, rotate, flip, etc.).

Administration

Contains the main image management functions (send image, delete image from archive, mark as findings-relevant, etc.).

Customizing the Toolbox

You can customize the toolbox according to your needs by clicking on . Then click on empty or filled toolbox fields to assign new functions.

By using the configuration bar that appears at the bottom of the toolbox, you can also add entirely new tabs to the toolbox, change their order, or assign certain modalities to the existing tabs.

3.3.2.6 Tool Region

The tool region allows you to store invisible functions in specific regions of the Viewer screen which can be quickly selected by clicking the middle mouse button (mouse wheel). To temporarily show the functions for orientation, press the **ALT + F1** shortcut and click into the Viewer. The tool region is available in all viewer windows.

Your administrator configures the available tools. The following default functions are available:

- Rotate 90° clockwise
- Scroll
- Rotate 90° counterclockwise
- Zoom
- Window/Level
- Pan

The mouse pointer only changes into the selected tool if said tool is also active in the current Viewer (e.g. scrolling is not possible in the multiframe Viewer).

The Tool region is available for all Viewers.

3.3.2.7 Presentation States (PR)

Presentation States (Grayscale Softcopy Presentation States) contain all the additional data relevant for image display, such as rotation, zoom, windowing, Lookup Table (LUT), graphical annotations, etc. These values are provided by the uploading modality and can only be partially changed in *syngo.share* view.

Opening a Study with Presentation State

Opening a DICOM study in *syngo.share* view automatically loads the associated Presentation States and applies the most recently created Presentation State. If multiple Presentation States exist, you can select the desired Presentation State by clicking the **PS** button.

If you have caching enabled, an exclamation mark may be displayed. The tool tip provides information about the error and how to fix it. These types of studies can be opened the same way as any other study. DICOM studies with associated Presentation States can be identified by the PR label in the Modality column.

Editing Presentation States


When loading a study with Presentation States to the Viewer, this study can be edited using applied Presentation States or non-applied presentation states (using windowing, zoom, pan, etc.). These changes can be undone by pressing the `Backspace` key to restore the original state. Changes cannot be saved.

3.3.3 General Functions

3.3.3.1 Assign Properties

☰ Copy view settings

Assigns the Zoom and Pan view properties to one or more images.


Note

If the `Viewing properties` `Assign view properties` function is enabled, the Zoom and Pan functions are applied to all images of the series.

☰ Copy level/window settings

Assigns the window leveling and pseudo colors view properties to one or more other images.

Note

If `Viewing properties` `Assign Window/Level` is enabled, the window leveling and pseudo colors functions are also applied to all images of the series. If the property is assigned afterwards, the mouse pointer changes and you can assign the property by clicking.

3.3.3.2 Calculations, Text Annotations, Drawing

syngo.share view lets you perform length, area, and angle calculations in images, use drawing tools for highlighting, and add textual annotations. The tools listed below are available for drawing and calculating areas, distances and angles, see [Section 3.3.3.2](#).

General Information

Shapes are drawn by holding the mouse button. Lines and polygons are drawn using single clicks for the starting and each intermediate point and a double-click for the end point. Arrows are drawn using single clicks for the starting and the end point.


Upon selecting a tool, the drawing toolbox opens. Select text size, line width and color from there.

A drawn object can be moved to any other image of a series within a layout by dragging it with the mouse pointer. Furthermore, a drawn object can be applied to all images of a series simultaneously by holding the `Alt` key while drawing. Please read the safety notices, see [Section 1.3](#).

Drawings can be transformed into calculations using the function `Show text for this annotation`, see [Section 3.3.3.2](#).

Calibrating Calculations

Manual calibration provides greater accuracy if a scale at level with the object to be calculated was included on image exposure.

**Note**

If calibration has not been performed yet, the displayed value refers to the calculated pixel spacing.

Please read the safety notices, see [Section 1.3](#).

Labeling Calculations

Calculations without calibration are labeled in pixels or square pixels. Calculations with calibration are displayed in millimeters or square millimeters. Depending on the use of pixel spacing, Presentation State or manual calibration, calculations receive additional labels:

det

Calculation is performed on the detector (Imager Pixel Spacing).

mag

Calculation is performed with magnification (Imager Pixel Spacing and Estimated Radiographic Magnification Factor).

ps

Calculation is performed based on the Presentation State.

man

Manual calibration has been performed.

us

Calculation is performed on an ultrasound image. See [Available Tools](#).

In all other cases, no additional labeling is performed.

Available Tools**3-point angle**

The 3-point calculation is intended for angles whose apex position can be determined in advance.

1. Place the crosshairs on a point of the ray forming the first leg and left-click.
2. Click a second time to set the apex of your angle.
3. Clicking a third time draws the second leg starting from the apex.

4-point angle


The 4-point calculation is intended for angles whose apex position is not known in advance.

1. Single-click a point on the first leg [1].
2. Click a second time on another point on the first leg. The straight line between points [1] and [2] is extended into infinity.
3. Next, click a third time to set a point at the second leg [3].
4. Possible apexes for the angle result from the point of intersection between the infinite straight line where points [1] and [2] are located, and the straight line where point [3] is located. The straight line with point [3] can still be moved and isn't fixed until you click a fourth time to set the apex [4].

Arrow


☰ Calibrate

To manually calibrate the tool using two points, click at the beginning and at the end of the scale in the image. Enter the length you drew according to the scale into the opened calibration window and apply. You can calibrate multiple segments by holding the **SHIFT** key. If a manual calibration was applied, all calculations are marked with (man).

Note

If a series is reloaded into the Viewer area, all calculations will be retained, but the calibration will be reset to the initial value. In addition, calibration changes are not saved when saving calculations.


Note

Anisotropic calibration for calculations on distorted images is not supported to prevent the user from applying this calibration to non-distorted images by mistake.

↳ Calibrating ultrasound images


Since ultrasound images show a curved area and may contain other image data besides the actual ultrasound, calculations are performed with the help of the ultrasound calibration region defined in DICOM. This region functions like a flat surface and can thus be used for calculations.

syngo.share view automatically shows the calibration region as a red frame when you select an annotation tool in the **View** area. Simply use the chosen annotation tool as usual to perform a calculation. Length and area calculations will be output in millimeters (mm) and square millimeters (mm²) and additionally labeled with (us).

Note

Please note that all points of your drawn line or shape must lie within the calibration region. Otherwise calculation is impossible and the calculation is output as "n/a". To correct this, move the drawing completely within the boundary set by the calibration region.

Note

It is theoretically possible for an ultrasound image to contain multiple calibration regions. In this case it is not possible to perform calculations via calibration region – you will have to calibrate manually instead.

Please read the safety notices, see [Section 1.3](#).

☰ Color

A submenu opens for the selection of color. You can either choose the desired color directly from the list or open the color selection window and pick a color from there.

☰ Ellipse

To draw a circle, press and hold the **Shift** key while drawing an ellipse.

☰ Freehand line

Hold the mouse button while drawing the line. To finish the line, release the mouse button.

☰ Freehand area

Press and hold the mouse button and trace the area of interest. (While drawing, the system always draws an interpolated line from the starting point over the supporting points back to the starting point.) If the lines intersect during or upon completion of the calculation, no area is calculated.

☰ Line/Length


A line can be a simple straight line or a continuous succession of multiple straight lines. The distance is calculated based on the most recent calibration. The calculated distance is shown next to the last point.

☰ Line width

Select the desired line width from the submenu.

☰ Lock annotation

This lets you permanently associate an annotation with an image (it is burned in). Therefore, you will be asked for confirmation before the action is performed.

Note

Locked annotations cannot be changed or removed as they are embedded in the document itself.

☰ Polygon

☰ Rectangle

To draw a square, press and hold the Shift key while drawing a rectangle.

☰ Remove all annotations

This removes all drawn objects and calculations. To delete a single drawn object, select it by clicking on one of its lines with the pointer. Press the Del key.

☰ Show text for this annotation


This function transforms a selected drawing into a calculation (and vice versa). The calculated values are shown above the drawn object.

☰ Text

Select the text size in the drawing tools window. Click on the area of the image where you want the text to appear, enter your text and apply it.

☰ Text size

Select the desired text size from the submenu.

Note

If you chose to calculate the drawn area, the calculated values are shown next to the shape. For CT images, density information is displayed as well:

- Average density (arithmetic mean value)
- Minimum density
- Maximum density
- Standard deviation in the density distribution

Editing a Drawn Object

To edit a drawn object, select it with the pointer by clicking on one of its lines. Change the length and/or shape of the object by dragging the white supporting squares which appear on the object. Dragging a drawing line (but outside the supporting points) will move the entire object without changing its size or shape. To move an object exactly horizontally or vertically, move it while holding the **Ctrl** key and the mouse button. You can also freely rotate a selected drawn object by using **Rotate and flip** **Rotate free**, see [Section 3.3.3.16](#).

When selecting a drawn object, the drawing tool opens in a separate window. Use the drawing tool to change:

- label text size,
- line width,
- line and labeling color,
- display of a label,
- permanent locking of a line or calculation.

You can perform and move calculations on multiple images of a series simultaneously.

Please read the safety notices, see [Section 1.3](#).

Deleting a Drawn Object


You can remove an annotation by clicking it with the pointer tool and then pressing the **Del** key.

Saving a Drawn Objects

Initially, drawings only exist on your PC. You have to save them together with the image. Saved drawings are then retained in the image or in a series and can be retroactively edited, but only if they have not been previously locked (burned into the image), see [Section 3.3.3.2](#).

As soon as a new annotation has been added, a **Save** function will appear in the related row of the Document Tree. This function will save the annotation to the archive.

When removing a document containing new drawings from the Document Tree, logging out or closing *syngo.share* view after drawings were done, the system will ask if you want to save the changes.

Note

Only annotations are saved in the archive. Calibration changes are not saved.

3.3.3.3 Camera

Transfers images, the complete Viewer screen or series within the Viewer to the print layout, see [Section 3.4.4](#). Please read the safety notices, see [Section 1.3](#).

☰ Photograph image

Point at the image you want to transfer to the print layout. Briefly press the left mouse button.

☰ Photograph Viewer

Point at any location in the Viewer to transfer the content of the Viewer to the print layout. Briefly press the left mouse button.

☰ Photograph series

Point at a series to transfer the complete series to the print layout. Briefly press the left mouse button.

3.3.3.4 Clear Viewer

Removes loaded elements from the Viewer. The respective elements remain loaded in the **Document Tree**. Either select **Clear Viewer** from the Viewer context menu, press the **Del** key or use the **X** button above the Viewer.


3.3.3.5 Density

☰ Show Density

Performs spot density calculations in CT images (e.g. bone density). Move the mouse across the region of interest to show the density in Hounsfield Units (*HU*) in a small field next to the cursor. Please read the safety notices, see [Section 1.3](#).

☰ Show average of density

Calculates the average tissue density in CT images for a circular region. The calculations use Hounsfield Units (*HU*). Click the center of the region whose average density you want to calculate and draw a circle while pressing and holding the mouse button. This circle marks the area for which the average density is calculated.

Note

This calculation is not saved. To save average density calculations, use the annotation tool.

Please read the safety notices, see [Section 1.3](#).

3.3.3.6 Image Administration

☰ Delete image from archive

Deletes the selected image from the archive (after a confirmation prompt).

☰ Send image

Sends the selected image to the archive. A window opens to select the target address (only DICOM-capable data recipients are supported). A status window shows the sending progress. Local images cannot be sent. For more information on sending data, see [Section 3.4.9](#).

3.3.3.7 Integration

This allows you to open external programs directly from *syngo.share* view. *syngo.share* view supports the integration of programs such as *mediCAD II*, *Voxar 3D*, or any other external applications. However, these external programs have to be configured in **Settings** > **Advanced** > **External Applications** beforehand.

Any changes made in an external application can be imported to the archive. Save the modified file in the respective external application, and then click **Save** in the Document Tree.

syngo.via and syngo.plaza

It is possible to send patient contexts to *syngo.via* and *syngo.plaza* in order to view documents quickly in either of these applications.

On startup, *syngo.share* view will check whether *syngo.via* or *syngo.plaza* has been installed and automatically configure the respective application in *syngo.share* view. External applications can also be manually configured under **Settings Area** > **External Applications**.

Text and Image Editing Software

Documents that cannot be displayed in the Viewer can be called up from the archive using an external application. In the case of a text file, for example, the appropriate text editing application will be opened.

Please read the safety notices, see [Section 1.3](#).

3.3.3.8 Invert

Inverts all grayscale values of an image (negative image).

3.3.3.9 Layout

This function is used for dividing the main window of the Viewer into multiple windows. Press and hold the mouse button and select as many fields as you want among the small squares displayed. Select only one field to use the main window for a single image. See [Section 3.3.2.2](#).

3.3.3.10 Mode Selection

Switch between the different modes available for the document type shown.

3.3.3.11 Pan

This tool () can be used to pan an image while pressing and holding the left mouse button.

3.3.3.12 Playback Controls

The playback control functions depend on the view selected. The functions available in movie mode (see [Section 3.3.5](#)) differ from those in the image series mode.

≡ First/Last Image

 and  show the first and last image of a series, respectively.

≡ Image selector

By pressing and holding the mouse button, you can use the slider to scroll through the series image by image.  and  serve the same purpose.

≡ Layout Step

The control buttons  and  are available when a layout with more than one section is selected in the display window in series mode. Next/previous layout takes you to the next set of images from that series, the step size depends on the layout chosen.

≡ Movie mode-specific functions

Detailed information on the various movie mode-only playback controls is available in [Section 3.3.5](#).

3.3.3.13 Pointer

This is the standard tool. It disables previously selected tools, selects any available calculation object and enables a tool in the tool region (see [Section 3.3.2.6](#)).

3.3.3.14 Print

Prints the displayed image, see [Section 3.4.4](#). Please read the safety notices, see [Section 1.3](#).

3.3.3.15 Pseudo Colors

This functions displays grayscale values in pseudo colors to highlight details. Click the desired option field to select a specific color scale. The first entry in the table is the default setting.

3.3.3.16 Rotate and Flip

Rotates or flips an image or a selected drawn object. The **Rotate and flip** submenu of the Viewer context menu contains the following tools:

☰ Rotate the image 90° clockwise or counterclockwise

☰ Rotate the image 180°

☰ Rotate free

This function rotates an image at a random angle. When you select the tool, the cursor changes to a circular arrow. Press and hold the left mouse button and freely rotate the image or use fixed angles:

- Turning the mouse wheel rotates the image in increments of 5°.
- **Ctrl + mouse wheel** rotates the image in increments of 45°.
- **Shift + mouse wheel** rotates the image in increments of 1°.
- **Ctrl + Shift + mouse wheel** rotates the image in increments of 9°.

☰ Flip

Flips the image along the horizontal or vertical axis.

If you flip an image horizontally, vertically, or both, the Viewer shows one of the following labels in the bottom right corner:

☰ S: h


for horizontal

☰ S: v

for vertical

☰ S: hv

for both


Note

Be careful when flipping images.

3.3.3.17 Scroll

This lets you scroll through the images of a series by pressing and holding the left mouse button and moving the mouse horizontally, or by turning the mouse wheel.

Note

If **Always use mouse-wheel for scrolling** is enabled in [Settings](#) > [General](#) > [Viewer](#), you can scroll using the mouse wheel without selecting the tool first.

3.3.3.18 Series Synchronization

☰ Synchronize series

Synchronizes the view in the Viewer windows.

☰ Synchronize series manually

This function connects two or more series in two or more Viewers based on the information of the loaded series (typically used e.g. for scrolling in matching series taken once with and once without a contrast agent). Scrolling in one Viewer entails synchronized scrolling in two or more Viewers.

This feature is particularly useful for slice images, such as CT scans. The slice thickness is taken into account during synchronized scrolling and displayed in the Viewer next to the ST (Slice Thickness) label.

The slice image series belong to the same study, the individual slice images are parallel to one another and have the same DICOM element `FrameOfReferenceUID` (can be checked using the **DICOM dump** function). The slice thickness (DICOM tag `SLiceThickness`) must be consistent within a single series, and no slice images must be missing. Different series, however, can have different slice thickness values.

3.3.3.19 Sorting Criteria

By default, documents are sorted by document number. This information is retrieved from the DICOM data and shown in the lower left corner of the image as the B: value. By default, the slice position is shown in the lower left corner of the image as the SL: value. A typical example for sorting a series by slice position is if the imaging device stores the document number incorrectly, causing the slice images to be displayed in the wrong order.

- Ascending sort (1, 2, 3, ...)
- Descending sort (... , 3, 2, 1)

The sorting order can be reversed by clicking the menu item.

3.3.3.20 Switch Screens

☰ Switch all Viewers to the next screen

With this function, all Viewers are moved over to the next screen in their current layout. This is especially useful when the second screen is located elsewhere (not directly next to the main screen) as all functions for viewing are available while the main screen remains free for other work.

☰ Switch all Viewers back to the previous screen

All Viewers are moved back to the previous screen in their current layout.

3.3.3.21 Tag

Tags the displayed image using either a keyword from a catalog or, if available, the keyword “findings-relevant”. Tagged images can be identified by the orange triangle. See [Section 3.2.4.4](#).

3.3.3.22 Viewing Properties

Always use interpolation (slow)

This function will recalculate an image with every editing step. This keeps the image quality high but requires more processing power. If this function is disabled, interpolation occurs when the image editing functions are completed.

Increase/Decrease text size

Incrementally increases/decreases the size of the text shown in the image (annotation text items are not affected).

Show annotation text

Determines whether annotation text items (e.g. from calculations) contained in the image are visible or not (the default text items of a DICOM file are not affected). Hidden annotation text items can be redisplayed at any time.

Show annotations

Determines whether calculations contained in the image are visible or not. Hidden calculations can be redisplayed at any time. To show/hide a single calculation, select it using the pointer tool first.

Show localizer

When scrolling multiple overlapping slice image series (CT and MR), you can display the position of the current slice in the other series. Marking the position of a slice in another series of slice images is called localizer display. The localizer display is disabled by default. To enable or disable the localizer display, select **Viewing properties** > **Show localizer** in the Viewer context menu. The term localizer is also used as a series description in scout images (also called topogram).

Depending on the type of series, the localizer appears slightly differently in the corresponding scout image. In the case of CTs and MRs, the localizer appears as a colored line – orange denotes the exact area shown in the series, blue denotes the level on which the area depicted is located.

Ophthalmic images, i.e. ophthalmic photography (OP) and ophthalmic tomography (OPT), are a special case due to the highly curved shape of the retina. For this reason, its localizer information is recorded in a special DICOM tag (Ophthalmic Frame Location Sequence). Depending on the type of series, the appearance of the localizer varies – if the ophthalmic image orientation is designated as linear, the localizer appears as an arrow; if it is designated as non-linear, the localizer appears as path segments. If multiframe series are hung, all corresponding localizers are displayed on the scout image and may overlap.

Show overlays

Enabling this option shows overlay pixel data contained in DICOM images.

Show text overlays

Determines if text items from a DICOM file shown on the image are visible or not (annotation text items are not affected). Hidden text items can be redisplayed at any time.

Synchronize Viewer settings


This function automatically applies certain viewing properties (e.g. zoom factor) for all images of a series in one Viewer (however, in case of a layout with multiple displayed items, in every single image). The following properties can be synchronized:

- Zoom
- Pan
- Windowing
- Auto level selection
- Pseudo colors

3.3.3.23 Undo Editing

Reset

This function undoes editing steps performed on the image (zoom, rotate, flip, etc.) and resets it. A series is reset to the first image.

Note

This does not remove inserted annotations (drawings, calculations, and comments).

Undo

Undoes the last action.

Restore

This function lets you restore the action that was just undone.

3.3.3.24 Level/Window (DICOM only)

The `Level/Window` function controls the brightness value of the image to increase or decrease the image contrast. You can change the brightness and contrast of the image display using the mouse.

Types of Windowing

The default windowing in the DICOM format is controlled using the Level/Center and Width settings.

syngo.share view also offers windowing using a lower and upper threshold (LT and UT). These thresholds restrict the brightness values displayed to the lower and upper limit of the respective histogram. Contrary to moving the width, LT and UT can be moved independently of each another. However, the values cannot overlap and can only be approximated up to their respective center value.

The preset window values for series which you open in the Viewer are always calculated across the entire series to ensure a uniform display. The window level is set as the mean value for all images contained in the series.

Lower and upper thresholds windowing is mostly used for images of nuclear medicine (e.g. NM, PT). Thus, *syngo.share* view can be configured to use an upper-/lower threshold windowing for these particular images by default (see `Settings` `Viewer`).

An option for non-linear windowing is VOI (Values Of Interests) LUT (Lookup table), see [Section 3.3.3.24](#).


Windowing Tools

☰ Level (brightness/darkness)

Press and hold the left mouse button and move the mouse horizontally to change the brightness.

☰ Width (contrast)

Press and hold the left mouse button and move the mouse vertically to change the contrast.

Note

With the center mouse button (mouse wheel), you can switch to the Brightness/ Contrast tool at any time without disabling the currently selected tool. Press and hold the center mouse button (mouse wheel) and move the mouse vertically and horizontally. When you release the center mouse button, the previous tool is once again selected.

☰ Advanced Level Tool ↪Auto Window/Level

Auto window/ level works independently of the recently used type of windowing (using Level/ Width or Lower/Upper Threshold) and shows all brightness values of the image.

↪Default value

This menu item resets the windowing of an image to the original value stored in the DICOM format. Other image editing steps are not undone (e.g. calculations, zoom, pan, etc.).

↪Skull, Lung, ..., Petrous bone

You can apply optimized presets for windowing images of specific body regions. Depending on the previously saved body region, these are available based on the device.

☰ VOI LUT


A VOI (Values Of Interests) LUT (Lookup Table) is a type of non-linear windowing.

If the DICOM images do not contain VOI LUTs, the VOI LUT button is not available. When loading DICOM images with one or more VOI LUTs into the Viewer, you can switch between the available VOI LUTs or display the DICOM images without VOI LUTs.

In addition, the windowing of DICOM images using VOI LUTs can be performed manually by compressing and stretching the VOI LUT curves.

☰ Auto level selection

After drawing a rectangular area in the Viewer, the system calculates the optimum windowing settings. The windowing settings are then applied to all images of a series. See [Section 3.3.3.24](#). This tool is only available for CT images.

Note

For DICOM images that do not utilize the maximum possible grayscale range, limits can be set in the *syngo.share* view configuration menu which restrict Level/ Width windowing. See [Section 2.2.4](#).

3.3.3.25 Zoom

☰ Zoom

Always keeps the focus on the point you click with the magnifier (precisely the center of the magnifying glass). To zoom in, press and hold the left mouse button and move the mouse up. To zoom out, press and hold the left mouse button and move the mouse down.

☰ Zoom selection

Lets you zoom in on a selected part of the image. Press and hold the left mouse button and draw a rectangle around the part you want to zoom in. Releasing the mouse button will zoom in the area of the rectangle.

☰ Magnifier


Enlarges a rectangular image section, which can be moved across the image while pressing and holding the mouse button.

☰ Fit to window

This function allows you to adjust the display of an image to the size of the current Viewer window. This is useful when you have changed the display for a single image (e.g. by enlarging it with a zoom tool) and want to quickly return to an optimal display size.

☰ Original size

This sets the zoom factor for the selected series to 100% and centers the images in the Viewer.

Note

Viewing properties > **Synchronize Viewer settings** is enabled by default, meaning that all images of the current series are zoomed 1:1, not just the selected image. You have to disable **Synchronize Viewer settings** to zoom only the currently selected image to 1:1.

The interpolation algorithm currently selected in **Settings** > **Viewer** will remain active.


3.3.4 ECG Mode

The ECG display is used to view waveform data from DICOM modalities.

In addition to several generic tools, the following special ECG functions are available:

☰ Caliper*

Adds visual interval guidelines. Click twice on the content to create an interval, which will be automatically applied across the data set.

Note

The visual lines created with the caliper tool will not be permanently stored in the document. They are used for visual purposes only.

☰ Horizontal Axis

Adjust the horizontal scale.

☰ Vertical Axis

Adjust the vertical scale.

☰ Data Layout

Arrange the way electrode data is displayed on the screen.

☰ Dataset Selector


Switch between datasets.

☰ Zoom In & Out (ECG Only)

Increase/decrease the zoom factor.

3.3.5 Movie Mode

Movie mode shows a series, e.g. multiframe images such as XA (X-Ray Angiography) or nuclear medical video sequences, as a film (multiframe view).

Note

When you load multiframe images into the series list, you do not have to wait until all frames are completely loaded. You can drag a series into the Viewer and view it while loading is still in progress.

The controls are similar to common music or video players, see [Section 3.3.3.12](#).

The context menu and the toolbox can be used the same way as in stack view.

However, the multiframe toolbar is specially designed for rendering multiframe images as movie sequences and is therefore radically different from the toolbar in stack view.

In addition to several generic tools, the following special movie functions are available:

☰ Play

 plays a movie sequence forwards and backwards in movie mode, and shows the next/previous image in series mode.

☰ Pause

 pauses the playback of the movie sequence. Use  to resume.

☰ Stop

 stops the playback of the movie sequence and jumps back to the beginning.

☰ Step

 and  the movie sequence forward/backwards by one step. The input field between these two buttons allows you to select how many individual images are included in one step.

☰ Replay Mode

The Replay Mode buttons provide different playback types for a movie sequence:

-  plays the movie sequence in an infinite loop.
-  plays the movie sequence in an infinite forward and backward loop.
-  plays the movie sequence forward one time.

☰ Scene Editing ↳Subscene

By clicking on , you set the starting point for a user-defined play sequence. Click again to set the end point of your sequence. Now the selected range is played repeatedly. Click the button again to end range playback.

↳Fix Scene

 is used to fix the scene, meaning that only the current scene is played repeatedly (the padlock is displayed closed). To end the current scene playback, click  again (the padlock icon is displayed open) so that all scenes in the sequence will be played back again.

☰ Multiframes with multiple scenes

Additional functions are available when multiframes contain multiple scenes:

↳ Scene quick-select

Clicking on allows you to quickly select and play a specific scene by clicking it (even during movie playback). If the movie has more scenes than shown in the window, you can use and to quickly switch to additional scenes.

↳ Movie scene selection

You can select a specific movie scene according to its number.

☰ Current frame rate


Use the input field to set the play speed. The set value indicates the number of frames per second. If you set a higher value than your computer can process, a warning symbol appears next to the frame rate.

If the movie starts automatically, the function **Autoplay multiframe images** is enabled in .

If movie mode does not start for multiframe series, stack view is selected for View multiframe series in under .

3.3.6 Video Mode

Video files (e.g. MPEG4) are automatically opened in video mode. The functions are similar to common video or music players, offering a rudimentary button and volume control.

Note

In case a video is not displayed properly, the video may have been imported with a variable bitrate.

In addition, you can add bookmarks to a video, meaning you can highlight and comment on particular sequences. By placing bookmarks, you can immediately jump to the desired part(s) of the video.

To bookmark:

1. Go to the desired part.
2. Click the button or press the key.
3. Enter your comment/description and click .

Using and , you can jump to/between different bookmarks. Alternatively, you can use the keys.

To delete a bookmark, use by clicking on the small arrow next to or press the key.

3.3.7 Metadata Mode

The metadata mode presents meta information of the selected element (e.g. DICOM image, Structured Report, etc.). You can use it to display and search the DICOM/EXIF Dump (i.e. DICOM or EXIF tags).

If a study contains more than one image, select the desired image from the selection list.

As long as the metadata view is active, every document that is loaded into the Viewer has its metadata displayed. To switch to other views (e.g. movie mode), use the button in the lower right corner of the Viewer toolbar.

3.3.7.1 Searching Metadata

You can search the metadata information for any text.

To find a specific text, enter the search term in the search field. As soon as two characters are typed in, the search will start immediately and automatically. The search is case-insensitive. If the search term is changed, the search automatically starts again. With the help of **F3**, more occurrences of the search term can be found. If you press **F3** again once all instances of the term have been found, you will receive a message stating that no further results could be found.

3.3.7.2 Copying Metadata

You can copy metadata to the clipboard and then paste it into another program.

Select the desired metadata information by pressing and holding the mouse button or select everything by pressing **CTRL + A**. You can also use **Select all** from the context menu.

Press **CTRL + C** or select **Copy** from the context menu.

The selected range of metadata information is now in the Microsoft Windows clipboard.

Open the program where you want to paste the information and press **CTRL + V** or select **Edit** **Paste** from the main menu of the program to paste the information from the clipboard.

3.4 Document Management

3.4.1 Merging Patients

Patient data may have to be corrected for a variety of reasons (e.g. because an emergency patient's identity was initially unknown or DICOM CDs were brought in from the outside and loaded into the system).

Patient corrections are typically performed by local administrators who are assigned to one or more organizational units. Using **Patient merge**, you can find inconsistent data and synchronize them using the in-house patient management system.

The tab **Patient merge** can be found in the **Admin** area in the main menu if you have the permission.

3.4.1.1 Finding Created Patients

In the **Patients** area, you can search for patients created by those organizational unit(s) for which you are the local administrator.

1. Select the timeframe in which you want to search for newly created patients.
2. Click **Refresh** to search within the selected timeframe.
3. The results list shows all patients created in the archive during said timeframe.

 allows you to export the results list as a table in HTML format to a text processing application of your choice or as a formatted table to Microsoft Excel.

All columns in the results list can be sorted by simply clicking the column header.

3.4.1.2 Finding a Matching Patient in the Patient Management System

In the search results list, click to select a patient whose information you want to verify.

This automatically starts a search in the patient management system for similar patients.

If search results were returned from the patient management system, they are listed in the **Auto matching patients** area.

If a link symbol is shown next to a result, it means that there exists a visit for this target patient whose visit ID is the patient ID of the source patient (i.e. the patient used to start the search).

Different colors indicate the type of (mis)match of the individual fields:

Green

The content of the result field exactly matches the content of the corresponding field of the selected patient. For an exact match, percentage points are assigned for the quality index.

Red

The content of the result field does not match the content of the corresponding field of the selected patient.

You can start a search for documents associated with a given patient by right-clicking and selecting **Search for documents of patient**. *syngo.share* view will switch to the archive search. If you want to delete a document (e.g. a study), choose the corresponding option from the archive search context menu.


If **Auto matching patients** does not list the desired results (or lists no results at all), you can start a manual search in the area of the same name.

All available search fields in the **Manual search** area are already filled based on the data of the selected patient. You can remove individual entries by clicking on the accompanying **X**, or, if you want to start a manual search from scratch, by clicking on **Clear**. **Search** starts the manual search according to your chosen criteria.

3.4.1.3 Determining the Correct Patient

To merge patients, click the correct patient in the search results list and then the **Set as Correct Patient** button (or the corresponding entry in the context menu).

For the time being, the system only makes note of the new link between the included documents and the correct patient. You can then repeat the patient correction for any number of patients.

Note

You can copy the content of a particular field in the search results list by right-clicking on it and selecting **Copy**.

3.4.1.4 Saving or Discarding Changes

Finally, you have to save all noted patient corrections:

1. Click the **Review and commit changes** button.
2. The list of new patients to be applied opens for verification purposes.
3. If all entries are correct, click **OK**.
4. If you find incorrect entries, click **Cancel**.

5. Then click .
6. Click at the confirmation prompt if you want to undo all the patient matches made.

3.4.2 Changing the Description

The following document properties and DICOM descriptions can be changed:

- [Study Level: Changing the Document Properties](#)
- [Series Level: Changing the Series Description](#)
- [Image Level: Changing the Image Description](#)

3.4.2.1 Study Level: Changing the Document Properties

In order to change a study or generic container description, examination date and/or time,

1. click the tab **Search**, or right-click on the study or generic container whose properties you want to change in the Document Tree.
2. Select from the context menu.
3. A dialog window opens where you can enter the new description, examination date and/or time.
4. Click .

3.4.2.2 Series Level: Changing the Series Description

In order to change a series description,

1. right-click on the series whose description you want to change in the Document Tree.
2. Select from the context menu.
3. A dialog window opens where you can enter the new description.
4. Click . The series description for the loaded series entry is updated.

3.4.2.3 Image Level: Changing the Image Description

In order to change an image description,

1. in the Viewer, right-click on the image whose description you would like to change.
2. Select from the context menu.
3. A dialog window opens where you can enter the new description.
4. Click . The image description is immediately updated in the Viewer.

3.4.3 Printing Documents

You can print various documents from the Viewer, including:

- PDF documents
- CDA documents
- Structured Reports (SRs)

Click on in the Viewer toolbar to open the print dialog. Select your desired printer, and click on . You can track the status of the print job in the **Job Queue**.

3.4.4 Printing Images

You can print DICOM and generic images from the Viewer. All calculation and layout settings are applied when printing images. The view currently set in *syngo.share* view is always used.

You can either print an image directly, or position/assemble one or more images in the **Print Layout** first.

3.4.4.1 Printing Directly

Click on  in the Viewer toolbar to open the print dialog. Alternatively, you can open the print dialog by selecting  from the Viewer context menu.

Select the desired printer, and make any other adjustments (e.g. changing the page orientation, adding a header or footer), and click on .


You can track the status of the print job in the **Job Queue**.

Changing the Background Color for Printouts

In the print dialog, you can change the background and text colors for the system printer and the DICOM printer. This way, you can delete black backgrounds and thus save on ink.


3.4.4.2 Printing via the Print Layout

In the **Print Layout** area, you can combine individual images, series and studies to print and export them as needed. Please read the safety notices, see [Section 1.3](#).


Note

If you add high volumes of data, it may take a few seconds before the result appears in the print layout.

Images in the print layout are always arranged across series by their image number. You can rearrange the images by clicking on and dragging them to the desired position. Once you are satisfied with the arrangement, click on  and follow the procedure described in [Section 3.4.4.1](#).

To remove images from the print layout, click on  and then click on the desired images. To resume arranging the remaining images, click on .

If you want to start compiling an entirely new print layout, you first have to click on . Otherwise, any new images you select will be added to the existing layout.

You have several options for adding images to the print layout:

- Dragging and dropping items from the Document Tree
- Photographing an image in the Viewer
- Photographing a series in the Viewer
- Photographing the complete Viewer

Dragging and Dropping Items

You can drag and drop the following items directly from the Document Tree to the print layout:


- single documents

- a series
- an examination
- all documents for a patient

To do this, proceed as follows:

1. Go to the **Print Layout** area.
2. Drag the desired series, examination, or patient from the Document Tree into the print layout by pressing and holding the left mouse button.

You can also select multiple specific items (series, examinations, patients) at once by pressing and holding **Ctrl**. For example, you can select two series of the same or different patients and add them to the print layout for a custom print compilation.

Note

If non-image files have been loaded for a selected patient, you cannot drag the patient to the print layout. Remove the non-image files from the Document Tree before proceeding.

Photographing an Image

To copy an individual image to the print layout, you can photograph it in the Viewer as follows:

1. Open the desired image in the Viewer.
2. Point the mouse pointer at the image you want to add in the Viewer.
3. Press **Spacebar** to photograph the image.
4. The photographed image is automatically added to the print layout.
5. Repeat this procedure for every image you want to add.

OR

1. Open the desired image in the Viewer.
2. Press **K** or select **Camera** > **Photograph image** from the Viewer context menu to activate the photography function (the mouse pointer changes).
3. Click on every image you want to add to the print layout.

Please read the safety notices, see [Section 1.3](#).

Photographing a Series

To copy the entire series to the print layout, you can photograph it in the Viewer as follows:

1. Open the desired series in the Viewer.
2. Point the mouse pointer at the Viewer.
3. Press **Ctrl + Spacebar** to photograph the series.
4. The entire series is automatically added to the print layout.

OR

1. Open the desired series in the Viewer.
2. Point the mouse pointer at the Viewer.

3. Press **Ctrl + K** or select **Camera** > **Photograph series** from the Viewer context menu (the mouse pointer changes).
4. Click on the Viewer to add the entire series to the print layout.

Please read the safety notices, see [Section 1.3](#).

Photographing the Viewer

You can also photograph the complete Viewer layout. Doing so adds all displayed images of a series to the print layout. For example, if four images (2x2 layout) are displayed in the Viewer, these four images are added to the print layout.

1. Open the desired series in the Viewer and select a layout other than 1x1.
2. Point the mouse pointer at the Viewer.
3. Press **Shift + Space** to add the displayed images to the print layout.

OR

1. Open the desired series in the Viewer and select a layout other than 1x1.
2. Press **Shift + K** or select **Camera** > **Photograph Viewer** from the Viewer context menu (the mouse pointer changes).
3. Click on the Viewer to add the displayed images to the print layout.

Please read the safety notices, see [Section 1.3](#).

3.4.5 Exporting Documents

Apart from opening data in an external application via the **Integration** function (see [Section 3.3.3.7](#)) `syngo.share` view allows to export documents in two more ways:

Export from the Print Layout

The data currently in the print layout tab are exported. To export data, click the **Export** button.

Export from the Document Tree

Select images, series, or all series from the Document Tree and select the export option from the document tree menu to export the data.

The export wizard opens. You can either export the selected data as DICOMDIR (see [Section 3.4.5.1](#)) or just save the file(s) in a selected format to a storage medium (see [Section 3.4.5.2](#)).

After starting an export process, you can check its status in the **Job Queue**.

Note

Threshold setting: If a threshold for loading the files is specified under **Settings** > **Advanced**, only this maximum amount of data is written when exporting. Similarly, setting a threshold can also be used to deliberately export only a portion of the series/container.

3.4.5.1 Exporting Documents as DICOMDIR


To create a DICOMDIR file, the wizard contains the following steps, depending on the options selected:

1. Export documents

- Select **As DICOMDIR** in the export wizard. You can now export a DICOMDIR to a hard drive. To burn the DICOMDIR to CD or export to a USB flash drive check additionally the appropriate checkbox to determine whether the file should be burned on CD.
- Check whether you want to anonymize patient data and/ or apply a DICOM export profile (see [Section 2.9](#)).
- Click .

2. DICOMDIR Settings

- Select the **Export with original transfer syntax** checkbox to prevent compressed studies from being decompressed during DICOMDIR export.

Note

Selecting this option may cause third-party Viewers not to be able to load the generated DICOMDIRs properly.

- Check the corresponding box if you wish to include *syngo.share viewLight* to the DICOMDIR.
- The File Set ID is an identifier for the generated DICOMDIR file designated by the DICOM standard. It may only contain uppercase letters, numbers, spaces, and underscores, and must contain 1-16 characters. If these conditions are not met, an exclamation point will be shown next to the field. This field is mandatory.
- The field for entering a directory is only available if you did not select the **Burn on CD** or **Write to USB medium** checkbox during the first step. Enter the drive and path where the file should be saved.

You can also click the button to select the directory from the folder selection dialog box or create a new folder if necessary. If you enter a path and the folder does not exist, it is created. Please follow Windows folder naming conventions (e.g. no special characters such as "\").

If available, this field is mandatory.

- Click .

3. Anonymization

- If you chose to anonymize patient data in the first step, you can freely change the specified fields or click the button to generate random data. Click .

4. DICOM export profile

- If you chose to apply a DICOM export profile in the first step, you can now select the profile and change the modification rules. Click .

5. Summary

- You see a summary of the actions you want to do. Click the button.

If you selected to save the DICOMDIR to your hard drive or to a USB flash drive, the export is completed. If you previously selected the **Burn on CD** option, the dialog box for configuring the CD burner opens.

6. Burning Settings

☰ Choose burn drive

If multiple CD burners are connected to your computer, you can select the desired device.

☐ with `syngo.share viewLight`

Select this checkbox if you want to burn `syngo.share viewLight` on the medium as well.

☐ Finish CD

It is recommended to select this option because multisession CDs which have not been finished cannot be read by all CD drives.

☐ Eject CD

Select this checkbox to eject the CD after burning is completed.

☐ ISO 9660 only

Selecting this option will burn according to the ISO 9660 standard, which means that the CD can be read on all operating systems.


After the burn process has started, the remaining and expired time are displayed alongside the burn rate.

3.4.5.2 Exporting Documents to a Storage Medium

To export data to a storage medium (without creating a DICOMDIR file) the wizard contains the following steps, depending on the options selected:

1. Export Documents

- Select the **Save to disc** option in the wizard.

Note

This method does NOT save the individual files as a directory but saves them individually in the selected folder.

- Click .

2. Export Settings

- Select the desired image format.

☐ Original file format (.dcm, .pdf, .avi)

Images are saved in the original file format and not decompressed.

☐ DICOM (.dcm)

Images are decompressed and resaved as DICOM.

☐ JPEG (.jpg)

Images are decompressed and resaved as JPEG. A selection list for the desired quality is displayed (High, Medium, Low). The higher the output quality, the more storage space is required for the exported files.

☐ JPEG2000 (.jp2)


Images are decompressed and resaved as JPEG2000. A selection list for the desired quality is displayed (High, Medium, Low). The higher the output quality, the more storage space is required for the exported files.

☐ PNG (.png)

Images are decompressed and resaved as PNG.

☰ TIF (.tif)

Images are decompressed and resaved as TIF.

Note

For documents that cannot be converted (such as PDF, txt...), only the original file format is available for export.

- Select the **Keep View settings** checkbox to export the images the way they are displayed in the Viewer.
- Select the **Create sub-folders for every single series** checkbox to save all images for a series in separate sub-folders. The sub-folders are called Series 1, Series 2, Series 3, etc., in the target directory.
- Select the desired numbering for the images from the selection list:

☰ Consecutive

Images are numbered in ascending order (0, 1, 2, 3, etc.).

☰ Document number/image number

Images are numbered according to the document/image number stored in the DICOM file.

☰ SOP Instance UID

Images are numbered according to the Service-Object Pair Universal Identification Number stored in the DICOM file. This is a unique 64-digit number that ensures that exported images receive a unique name.

☰ None

Images are not numbered.

- Select the **Front** checkbox to add the numbering as a prefix to the file name. If the checkbox is not selected, the numbering is added as a suffix to the file name.
- **Base filename:** This is the future filename, without numbering or file extension. It remains the same for every image of an exported series. You can choose from a selection list of standard combinations of DICOM tags or create your own names and add them to the selection list.

☰

Opens a new window where you can create your own base filenames to be added to the selection list for future quick retrieval. Click to save the new base filename.

☰

Deletes the selected base filename from the selection list.

☰

Opens a new window where you can edit the selected base filename.

☰

Resets the base filename in the field and the selection list to the default values.

- **Export Folder:** State the target folder or open the folder dialog box to select the target folder for the data to be exported.

- Click .

3. Advanced Settings: If you have chosen JPEG, JPEG200, PNG or TIF as the file format, you now have further options to customize the export:

- You can change the background and text color for the images to be exported. Enable or disable the various checkboxes to determine which image elements should be included in the export. The preview in the dialog box is updated according to your settings.
- Selecting the **Automated alignment** checkbox will automatically adjust the printout with regard to the labeling. If the checkbox is not selected, you can set the resolution in dpi, the image height and width (except for DICOM images).
- **1:1** exports the images in their original size. This option is only available when not exporting in DICOM format. Click on .

These settings apply to all images of the current export job. Click .

4. Summary

- A summary of your settings is displayed. Click the button to perform the action as listed in the summary.

3.4.5.3 Errors

⚠ DICOM directory available

If you selected a directory which already contains DICOM directory data, an error message is displayed.

You can append the data, overwrite them, or cancel the export.

⚠ Invalid DICOM tag value

Sometimes during export, DICOMDIR may include files containing an invalid value. This results in an error message.

You can display a detailed error message which may indicate where exactly the error occurred. Alternatively, you can select to add the file regardless or select to not add the file.

⚠ Not all documents were exported

You may encounter a message stating that not all documents were exported because generic documents (e.g. PDFs) could not be converted into DICOM format.

Change the export settings or the selected files if necessary and repeat the export procedure.

⚠ No documents were exported

This error message warns that no documents were exported and that the temporary DICOMDIR file will be deleted.

If necessary, change the export settings or the selected files and repeat the export procedure.

3.4.5.4 Creating a Dossier


The **Dossier** function lets you collect and export documents of one or multiple patients together with their relevant metadata and forward them to a recipient.

The selected files are automatically collected and placed on a preconfigured network drive in an encrypted ZIP file. This ZIP file contains a separate folder for each patient. Each patient folder in turn con-

tains a folder for DICOM and Generic Documents, as well as a CSV file with additional information on the patient (e.g. social security number, name, patient ID) and the included documents.

You can create a dossier in the **Archive** tab. To add a document, select it in the search results list, right-click and go to **Dossier** > **Add**. You can view and edit the dossier anytime via **Dossier** > **Show...**. The column **Already used** indicates if a document was previously added to a different dossier.

To create/forward the dossier, enter the recipient and a meaningful description in the relevant fields. You also have to define a password in order to protect the ZIP file against unauthorized access. Then click on **Create** to finalize the dossier.

Note

You can only work on one dossier at a time. If you want to create a new dossier, you have to first forward the current one.

You can check on the status of the dossier in the **Job Queue** (☰ in the main menu). Since some of the collected files may have to be retrieved from long-term archives or other VNAs, the creation of the dossier may take quite some time. However, you can continue working with *syngo.share* view as usual.

Once the dossier is finished, the recipient will receive an e-mail containing the name and description of the ZIP file. Since the file will only be available on the network drive for a limited period of time, you may want to download the dossier as soon as possible.

3.4.6 Anonymizing Documents

3.4.6.1 Anonymizing Documents during Export

Depending on permissions, documents may have to be anonymized during export.

When anonymizing, patient-relevant DICOM fields are filled with anonymized values. The content of optional patient-relevant fields is also removed (in 0010 DICOM tags). The following list shows the fields that are overwritten during the anonymized export (e.g. with automatically generated random values):

- PatientID
- PatientName
- PatientSex (optional)
- Patient Birth Date (optional)

If you anonymize the fields yourself or generate random values, or only have permission for anonymized export, the content of the DICOM OriginalAttributesSequence in the DICOM file as well as private tags are removed. Additionally, the filename that is created during export is anonymized.

3.4.6.2 Anonymizing Documents during Printing

It may be necessary to anonymize data when printing images. You have the following options:

☰ Show Text Overlays

You can hide image data in the viewer by pressing **T** or by clicking on the corresponding icon in the toolbar. The data stays hidden when printing.

Hide the patient data

By clicking on the icon in the main menu, patient data such as name, sex, etc. will be replaced by asterisks. This is a global function, i.e. even if you use **Show text overlays** and/or enter an image caption in the print menu, the data remains anonymized during printing.

Hide Image Caption in the Print Dialog Window


You can select **None** from the dropdown menu in the image description section of the print dialog window. However, if you select **As camera** instead, this will copy the image caption as shown in the preview. For more information, see [Section 2.8](#).

3.4.7 Moving Documents

Moving documents means changing the producer, i.e. the organizational unit, of an examination (DICOM study or generic document). This is a special feature in *syngo.share* view.

To move a document, proceed as follows:


1. Select a study in the search results list.
2. Right-click to open the context menu, and select **Move...**.
3. Doing so opens the list of possible recipients.
 - You can limit the displayed recipients by using the filter function.
 - If you have added recipients to your list of favorites, they will be shown separately (see [Section 3.4.15](#)).
4. Click to select the desired recipient from the list. To select multiple recipients at once, press **Ctrl** or **Shift** while clicking.
5. Confirm your recipient(s) with **Select**. Doing so moves the document.

Note

The results list is not updated until you run a new search.

After starting the process, you can check its status in the **Job Queue**.

Note

If a recipient is not listed, please contact your system administrator. This action can be carried out simultaneously on several documents as well. Simply select the desired elements and then perform the action.




See [Section 3.4.15](#) for managing quick access to recipients.

3.4.8 Copying Documents


Documents can be copied to another storage area (=organizational unit). The command is available in the context menu of the results list.

All calculation and layout settings are applied when copying images. The view currently set in *syngo.share* view is always used. Copying will insert the images in their original size. The size can later be adjusted in the target application.

To copy a document:

1. Select the document to be copied from the search results list.
2. Right-click to open the context menu, and select .
3. Doing so opens the list of possible recipients.
 - You can limit the displayed recipients by using the filter function.
 - If you have added recipients to your list of favorites, they will be shown separately (see [Section 3.4.15](#)).
4. Click to select the desired recipient from the list. To select multiple recipients at once, press  or  while clicking.
5. Confirm your recipient(s) with . Doing so copies the document.

After starting the process, you can check its status in the **Job Queue**.


Note

If a particular recipient is not in the list, please contact your system administrator. The process can be executed for several documents at once. Simply select the desired elements and choose the command.

3.4.9 Sending Documents

All DICOM-capable applications or devices integrated with the system can receive documents. Select the recipient from a list of *Application Entity Titles (AET)*. The AET is an identifier assigned to the receiving device by the system administrator in *syngo.share* Configuration.

You can send generic documents, studies, series and individual images.

Note



If a recipient is missing from the list, please contact your system administrator. The process can be executed for several documents at once. Simply select the desired elements and call the command. A simultaneous sending of DICOM files and generic documents is not possible. The recipients have to be configured differently.

See [Section 3.4.15](#) for managing quick access to recipients.

After starting the process, you can check its status in the **Job Queue**.


Sending DICOM Studies, Series or Generic Documents

To send generic documents, studies or series:

1. Select the desired document or study in the search results list or in the Document Tree. A series can only be selected in the Document Tree.
2. Right-click to open the context menu, and select  .
3. Doing so opens the list of possible recipients.
 - You can limit the displayed recipients by using the filter function.
 - If you have added recipients to your list of favorites, they will be shown separately (see [Section 3.4.15](#)).
4. Click to select the desired recipient from the list. To select multiple recipients at once, press  or  while clicking.

5. Confirm your recipient(s) with **Select**. Doing so sends the document.

Sending an Image

Note

You can only send one image at a time.

To send an individual image:

1. Open the desired image (DICOM or generic) in the Viewer.
2. Right-click on the displayed image to open the context menu, and select **Image administration** **Send image...**.
3. Doing so opens the list of possible recipients.
 - In either case, a list of possible recipients opens. You can limit the displayed recipients by using the filter function.
 - If you have added recipients to your list of favorites, they will be shown separately (see [Section 3.4.15](#)).
4. Click to select the desired recipient from the list. To select multiple recipients at once, press **Ctrl** or **Shift** while clicking.
5. Confirm your recipient(s) with **Select**. Doing so sends the image.

3.4.10 Sharing Documents

By sharing documents, you make them available to private practice physicians and medical institutions for electronic report communication.

To share a document:

1. Select the study or generic document you want to share in the search results list.
2. Right-click to open the context menu, and select **Share...**.
3. Click on **To** to select individual recipients or on **To groups** to select predefined groups of recipients. You can also combine individual recipients and groups.
 - In either case, a list of possible recipients opens. You can limit the displayed recipients by using the filter function.
 - If you have added recipients to your list of favorites, they will be shown separately (see [Section 3.4.15](#)).
4. Click to select the desired recipient from the list. To select multiple recipients at once, press **Ctrl** or **Shift** while clicking.
5. Confirm your recipient(s) with **Select**.
6. Select how long the share should be valid. If you do not select a validity period, the share will be valid for a default period.
7. Enter a subject and message for the recipient.
8. Optionally, you can assign priority to the share message by checking the **High Priority** checkbox.
9. Click on **Send** to share the document.

3.4.11 Registering Reference Pointers

You can manually register documents in order to make them available to subsystems. This means that a reference pointer containing all relevant information about the document(s) is sent to the subsystem.

To register a study or generic document:

1. Select the desired element in the search results list or in the Document Tree.
2. Right-click to open the context menu, and select **Register reference pointer...**.
3. A window with a list of possible recipients opens.
 - You can limit the displayed recipients by using the filter function.
 - If you have added recipients to your list of favorites, they will be shown separately (see [Section 3.4.15](#)).
4. Click to select the desired recipient from the list. To select multiple recipients at once, press **Ctrl** or **Shift** while clicking.
5. Confirm your recipient(s) with **Select**. Doing so registers the reference pointer for the study/ generic document.

3.4.12 Reassigning Documents

By reassigning documents, you can create the correct patient context for incorrectly saved documents.

Reassigning a document means relinking a document:

- To a patient with a visit
- To a patient without a visit
- To another visit
- From one entity to another
- From an entity to a patient and vice versa
- From an entity to a visit of a patient

To change a patient, a visit or an entity for a document:

1. Right-click on the element you want to reassign in the search results list. Select **Reassign to another patient/visit...** or **Reassign to another entity...** from the context menu.
2. A window opens showing searches for visits, and patients, or – in a separate tab – for entities.
3. Enter either the visit ID (AN), the data for a patient, or the data of the entity to whom/which you want to reassign the document.
4. Click **Search** to search for the visit ID, the patient or entity.
5. Depending on the search criteria, the result returns a list of one or all patients with visits or entities matching your search criteria.
6. If one or more visits exist for a patient, you can view them by clicking the **+** symbol.
7. Select the patient or a specific visit to which you want to reassign the examination.
8. Click **OK** to start the process.

3.4.13 Deleting Documents

The following documents can be deleted from the archive:

■ DICOM study or generic document

in the results list of the search or in the Document Tree

■ DICOM series


in the Document Tree

■ DICOM image or generic image

in the Viewer

When deleting a document, a dialog box opens where you have to select a reason for the deletion from the drop-down menu. The available reasons for deletion depend on the settings made by your system administrator.

You can delete several documents at once. In this case, the mandatory reason for deletion applies to all selected files. To delete multiple files, select them with **Ctrl + click** or **Shift + click** and execute the command. Note that you cannot simultaneously delete DICOM files and generic documents in the Document Tree.

Note

Caution: You cannot restore deleted documents in *syngo.share* view. This is only possible for your system administrator.

After starting the process, you can check its status in the **Job Queue**.

Deleting Generic Documents, Studies or Series

This is how you delete a generic document, a study or series from the archive:


1. Select the element you want to delete in the results list or in the Document Tree. A series can only be selected in the Document Tree.
2. Right-click to open the context menu and select **Delete**.
3. A security pop-up appears. Select the reason for deletion from the drop-down menu.
4. Click **OK** to delete the element.

The current results list and the Document Tree are updated immediately.

Deleting an Image

This is how you delete an image from the archive:

1. In the Viewer, right-click on the image you want to delete.
2. Select **Image Administration** > **Delete Image from Archive**.
3. A security pop-up appears. Select the reason for deletion from the drop-down menu.
4. Click **OK** if you really want to delete the image.

**Note**

Only one image can be selected at a time.

The deleted image is removed from the Viewer and deleted in the archive.

3.4.14 Viewing and Correcting Warnings and Errors

syngo.share view provides a document error and warning log . An error log is created when an error occurred during the import of DICOM files. In case corrective measures set in the DICOM configuration are taken during a re-import, a warning log is created.

To find files with an error or warning, use the search filter **Document Log**.

Follow these steps to access the document log and correct entries:

1. Right-click on a document in the results list to open the context menu and select **Show document log** .
2. A new window opens and lists all warnings and errors logged for the selected file. Warnings are marked with a yellow exclamation mark, errors are marked with a red exclamation mark.
3. Click on the desired entry to see the entire log entry.
4. If the warning is acknowledged or the error is resolved, e.g. after a re-import, you can delete it from the list. Select it and click **Set as resolved** . You can select several entries and resolve them all in one step. For traceability you may also keep the entries in the document log.

To resolve issues use the command line tools available. For more information please see the respective chapter in the OnlineHelp *syngo.share* core.

3.4.15 List of Favorites

The list of favorites usually contains preconfigured entries of preferred recipients (it may initially be empty). You can configure the list according to your needs. Click on the corresponding star symbol to add or remove a recipient to/from your list of favorites.

Alternatively, you can add/remove recipients via the context menu (right-click and select **Add to favorites** or **Remove from favorites**).

A list of favorites is available for the following actions:

- [Sending Documents](#)
- [Moving Documents](#)
- [Copying Documents](#)
- [Sharing Documents](#)
- [Registering Reference Pointers](#)

3.5 Keyboard Shortcuts *syngo.share* view

Search

F5	Start search
F2	Change document properties
F8	Toggle series thumbnails
F9	Load selection to document tree
Ctrl C	Copy selected search results (in HTML format)
Del	Delete document/generic container

Document Tree

F2	Change document properties
F4	Clear document tree
F7	Toggle document tree
F12	Show next patient in document tree
Shift F12	Show previous patient in document tree
Home	First entry in document tree
End	Last entry in document tree
Ctrl S	Export selected series
Ctrl Shift S	Quick export selected series

Viewer: General

F11	Toggle fullscreen
Alt F11	Toggle white screen
Shift F11	Toggle multiscreen
Tab	Move viewer grid to next screen
Shift Tab	Move viewer grid to previous screen
Ctrl Shift B	Change description (generic files)
Ctrl Shift D	Delete image from archive

Viewer: Navigation and Actions

→	Next image in series
←	Previous image in series
Home	First image in series
End	Last image in series
Ctrl →	Next viewer
Ctrl ←	Previous viewer
Ctrl ↓	Viewer down
Ctrl ↑	Viewer up
Page ↓	Next set of images (Layout view)
Page ↑	Previous set of images (Layout view)
S	Synchronize viewers manually
Alt S	Deactivate synchronization
Shift S	Synchronize all viewers
Y	Synchronize viewer settings
Shift Y	Apply viewer settings to all viewers
Del	Remove image from viewer
Ctrl Shift L	Clear all viewers
Ctrl C	Copy image

Viewer: Print and Export

Ctrl P	Print
Space	Load image at cursor position directly to print layout
Shift Space	Load entire viewer directly to print layout
Ctrl Space	Load entire series directly to print layout
K	Camera tool: Load active image to print layout
Shift K	Camera tool: Load entire viewer to print layout
Ctrl K	Camera tool: Load entire series to print layout
2-6	Open file in external application 2 to 6
Ctrl Shift T	Send image

Viewer: Visibility of Information

Ctrl L	Toggle localizer
O	Toggle overlays
T	Toggle labels
Alt T	Toggle annotation texts
Ctrl T	Toggle annotations
Ctrl Shift A	Show/hide patient data in viewer

Viewer: Tools

Esc	Pointer tool
B	Scroll through series
D	Tissue density calculation
F	Window
Ctrl Alt A	Auto-window
Alt F	Auto-window based on selection
Shift F	Default windowing value
Shift 1-9	Windowing presets
E	Ellipse
P	Polygon
R	Rectangle
L	Length
V	Pan
Z	Zoom
M	Zoom selected area
1	Zoom tool: Original size
Shift Z	Zoom tool: Fit to window
Ctrl +	Increase label size
Ctrl -	Decrease label size
Ctrl Z	Undo
Ctrl Y	Redo
Backspace	Reset

Video Viewer

Space	Play/Pause
→	1 second forward
←	1 second backward
↓	30 seconds forward
↑	30 seconds backward

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ITH icoserve technology for healthcare GmbH
Siemens Healthcare GmbH