syngo.share / Release VA30C / 2021-12-16 / Revision 6272

syngo.share webview

User Manual

Language Directory

Dansk	3
Deutsch	13
Ελληνικά	23
English	33
Español	43
Français	53
Italiano	63
Magyar	73
Nederlands	83
Norsk	93
Polski	103
Português (Brasil)	113
Português (Portugal)	123
Suomi	133
Türkçe	143
简体中文	153

syngo.share / Release VA30C / 2021-12-16 / Revision 6272

syngo.share webview

Brugermanual

Indhold

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ITH icoserve technology for healthcare GmbH
Innrain 98
A-6020 Innsbruck
Austria

2021

Distributed by

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens-healthineers.com

Legal Manufacturer

ITH icoserve technology for healthcare GmbH
Innrain 98
6020 Innsbruck
Austria
Phone: +43 512 89059

Made in Austria

syngo.share / Release VA30C / 2021-12-16 / Revision 6272

syngo.share webview

User Manual

Contents

1	Introduction	35
1.1	Intended Use	35
1.2	Hardware Requirements	35
1.3	Software Requirements	35
1.4	Additional System Requirements	36
1.5	Safety Notices	37
	Image and Result Distribution	37
	Medical Education	37
	Accuracy of Calculations	37
	Units	37
	Access Restriction	37
	Effects of External Programs	37
	Printing of Images	38
	Network Connection	38
	Workstation Accessibility	38
2	Usage Instructions	39
2.1	Access the Application and Log in with Your User Credentials	39
2.2	Search for Patients and Documents	39
2.3	Customize the Results List and Preview the Results	39
2.4	Manage Documents and Metadata	40
2.5	View Documents	40
2.6	Share Documents via the Inbox	40
2.7	Archive Documents	41

1 Introduction

1.1 Intended Use

syngo.share webview is a web application to be used in clinical image and result distribution and provides access to clinical data from connected systems (such as Vendor Neutral Archives, Picture Archiving and Communication Systems, DICOM archives, XDS-I Affinity Domains, Electronic Medical Record Systems, Clinical Information Systems) and the local file system for the electronic display of data. This product is not intended for diagnosis, prevention, monitoring, prediction, prognosis, treatment, cure or alleviation of diseases or other medical conditions.

1.2 Hardware Requirements

Table 1: Hardware Requirements syngo.share webview Client

Part	Requirements
Processor	≥ 2 Cores
Memory	≥ 4 GB
Storage	≥ 80 GB
Network Interface	≥ 100 Mbit/s
Network Delay	≤ 100 ms; ≤ 60 ms (recommended)
Recommended Screen Resolution	≥ 1280x1024

1.3 Software Requirements

- Windows 10 Enterprise / Professional Edition x64
- Windows 10 Enterprise / Professional Edition x86

The following web browsers are supported with an optimum resolution of ≥ 1280x1024:

- Google Chrome ≥ 83
- Microsoft Internet Explorer 11 (not recommended)
- Microsoft Edge ≥ 83
- Mozilla Firefox ≥ 78
- Mozilla Firefox ESR ≥ 78
- Apple Safari ≥ 13

syngo.share webview has been tested with the following web browsers:

- Google Chrome 89 (recommended)
- Microsoft Edge 89

Additional requirements for the backend:

- Red Hat Enterprise Linux 8 (≥ 8.2) (recommended) / Red Hat Enterprise Linux 7 (≥ 7.4)

Mobile Devices

The following mobile web browser is supported:

- Apple Safari (iPadOS) = 13, 14

syngo.share webview has been tested with the following mobile device and web browser:

- Apple iPad (2020)
- Apple Safari (iPadOS) 14 (recommended)

Feature restrictions may apply on mobile devices (e.g. export, download and print). Mobile view is required, touch input is recommended. It is further recommended to launch *syngo.share* webview via home screen shortcut ("website icon"). This will cause the browser to open in full-screen mode which offers the best user experience of *syngo.share* webview. (How to set a shortcut on the home screen is described in the online iPad user guide.)

1.4 Additional System Requirements

Note Regarding Videos

To allow for the highest possible loading speed when viewing videos, the following is recommended:

- The HTTP servers should be activated on FilesystemServers where video data is stored.
- Videos should be imported to an organizational unit for which an archiving rule without data compression is defined.
- The video codec H.264 (Level: 3.0, Profile: Baseline or Constrained Baseline) and audio codec AAC should be used.

Note Regarding Whole Slide Images for Digital Pathology

DICOM Whole Slide Images for Digital Pathology are only supported in conjunction with *syngo.share* core and have to be imported to an organizational unit for which an archiving rule without data compression is defined.

1.5 Safety Notices


Image and Result Distribution

syngo.share webview is only intended for clinical image and result distribution.
Do not use *syngo.share* webview for direct diagnoses.

Medical Education

The use of *syngo.share* webview should be granted to qualified staff members only.

Accuracy of Calculations

Images may contain calibrations. Please note that the accuracy of these calibrations correlate to the precision of the program or modality with which they were created. Calibrations must be checked according to the manufacturer's instructions. Calculations, calibrations and results derived therefrom must be checked for plausibility by clinical professionals. Inaccuracies may arise during image capture by clinical professionals or due to the image capturing technology of a given modality. For instance, X-ray images may contain inaccuracies due to the fact that the x-rayed body part is projected onto a two-dimensional imaging surface. Inaccuracies during visualization may be due to browser or hardware settings.

This product rounds calculations as follows: Length and area: 1 decimal place; density (HU): no decimal place; area density (HU): 1 decimal place. This product uses an interpolation algorithm for image display.


Units

The application *syngo.share* webview conforms to Directive 80/181/EEC, which contains specifications related to units of calculation. Calibrated lengths are calculated using millimeters. Surface areas are calculated in square millimeters (mm²). Plane angles are calculated in degrees. Densities are calculated using the industry-standard unit, Hounsfield units (HU), which are not regulated by Directive 80/181/EEC.

Access Restriction

Be sure to log out of *syngo.share* webview when finished working in order to avoid unauthorized access to the system.

Effects of External Programs

The software manufacturer does not guarantee that third-party applications will not negatively affect the system. The standards of any external programs accessed by *syngo.share* webview are set by their respective manufacturers.

Printing of Images

To print *syngo.share* webview transfers the data to the print function of the browser. Thus there may be differences between the view in the Image Display Area in *syngo.share* webview and the printout, e.g. the image could be cut off, annotations could be incomplete or missing entirely. Possible transformations such as mirroring, rotation, and applied presentation states are not marked explicitly in the printout.

Network Connection

syngo.share webview requires an active network connection. For IT security reasons, a (temporary) failure of the network connection may result in you being logged out and having to log in again.

Workstation Accessibility

When using *syngo.share* webview, make sure that you can operate the application at all times.

2 Usage Instructions

syngo.share webview is a web application that allows you to view and distribute medical documents (e.g. images, reports, videos and photographs). It offers a variety of built-in viewing and annotation tools. You can upload documents directly to a patient context, as well as share and send them. In addition, you can edit and manage patient and document metadata.


2.1 Access the Application and Log in with Your User Credentials

Provided you have not been logged in by a third-party application, use the following URL template to access *syngo.share* webview:

`https://<server>/webview/<tenant>/<site>`


Enter your user credentials to log in to *syngo.share* webview. If applicable, you also have to select the domain. If you want to preselect the domain, you can provide this information in the URL according to the following template:

`https://<server>/webview/<tenant>/<site>/?domain=<domain>`

Note

If your institution only has a single tenant and/or site, you do not have to provide `<tenant>` and/or `<site>`. If you are using *syngo.share* webview without *syngo.share* core, tenant and site do not apply.

Note

In case trusted-only login is activated for your institution, access to *syngo.share* webview is only possible via a valid SAML token or HMAC signature.

2.2 Search for Patients and Documents

To begin a search, enter the desired search criteria in the appropriate fields and press the button. The button will remain inactive if the entered search criteria are insufficient.

You can customize the visibility of search fields according to your needs and save search queries with predefined criteria as filters. Moreover, you can favorite a filter, which will execute it automatically upon login.

You can search for patients and documents in *syngo.share* core as well as other available DICOM nodes. You can also start a cross-enterprise search (XDS search) for and load external documents for a given patient, provided that an eHealth Solutions Consumer Adaptor has been configured for this purpose.

2.3 Customize the Results List and Preview the Results

You can arrange the various columns as desired and change the sorting order of the results. You can also group the results according to patients and documents as well as filter them via text entry.

Your adjustments are saved in your user profile and will be applied to all future search queries.

Selecting an entry from the results list provides you with detailed series information and referenced documents as well as an interactive preview. Various display options are available.

In addition, you can load available prior examinations for found DICOM document entries and preview their contents. Opening a prior examination in the viewer displays it side by side with the corresponding current examination.

2.4 Manage Documents and Metadata

Depending on your user permissions, various functions are available for managing documents, patients and their metadata. These functions can usually also be accessed on the series and instance level:

- Open documents in an external application
- Add, edit or delete metadata and keywords
- Export documents (anonymized export is available)
- Send documents to remote peers (e.g. DICOM node)
- Share documents to other applications (e.g. Patient Portal, Inbox, EMR systems)
- Copy documents, move them to another destination or assign them to another patient/visit
- Mark documents as “deleted”
- Rearrange studies and series

2.5 View Documents

Double-click on a document to open and display it in a predefined layout. You can customize the layout according to your needs.

Different viewing and annotation tools are available depending on the file type:

- Panning, windowing and scrolling functions
- Various transformation and zoom tools
- Various length, area and angle annotation tools
- Drawing and text annotation tools
- Image presets and display Presentation States and DICOM Overlays
- Add, edit or delete metadata and keywords
- Playback tools in *Cine Mode*
- Synchronization of series
- Display of Reference Pointers

The patient banner is always visible in the Viewer to prevent patient mix-ups.

The *Patient Jacket* displays all available documents in the archive for the current patient. Corresponding documents from third-party archives can also be loaded.

2.6 Share Documents via the Inbox

The inbox is used for the electronic transmission of documents and is available to all registered doctors.

You can share documents with one or more selected users. The recipient(s) are automatically informed when documents are shared with them.

The inbox offers the same functions as the results list (see [Section 2.3](#)) and the Viewer (see [Section 2.5](#)).

2.7 Archive Documents

syngo.share webview allows you to archive documents directly to patients or visits. Archiving documents is only supported for *syngo.share* core.

To archive documents, follow these steps:

1. Set a patient context.
2. Select the required data.
3. Select an archiving destination.
4. Add keywords.
5. Complete archiving.

ITH icoserve technology for healthcare GmbH
Innrain 98
A-6020 Innsbruck
Austria

2021

Distributed by

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens-healthineers.com

Legal Manufacturer

ITH icoserve technology for healthcare GmbH
Innrain 98
6020 Innsbruck
Austria
Phone: +43 512 89059

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